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## Q-LOGIC 2 BLUETOOTH OWNER'S MANUAL UPDATE

To add a bluetooth screen (Aux menu must be enabled by provider to use bluetooth functions):

1. Push the On/Off and Mode Select Lever up once to power on the chair and controller.
2. Push the Menu button to access the Main Menu screen. **See figure 1.**
3. Use the joystick to highlight the Handcontrol settings option, then push the joystick right to select that function. **See figure 1.**
4. Use the joystick to highlight the Bluetooth option, then push the joystick to the right to select that function.
5. Use the joystick to highlight the bluetooth function you wish to enable, then push the joystick to the right to select that function.

**NOTE: The first time you select each Bluetooth function, the "Powercycle Needed" message will appear. To proceed, you must use the On/Off and Mode Select Lever to power off, then power back on the controller to enable your selected Bluetooth function.**

6. Push the joystick forward to enable the bluetooth screen. **See figure 1.**
7. Push the joystick left to select the esc option and return to the bluetooth menu.
8. Press the menu button to return to drive mode.


Figure 1. Bluetooth Menu

**NOTE: This product has been tested and has passed Electromagnetic and Radio Frequency Inference. (EMI/RFI). Please refer to the Consumer Safety Guide for more information regarding EMI/RFI. Although we do not recommend using a cell phone while driving a power chair, the system is capable of receiving and/or making calls through a cell phone device, accessing music files, contact lists, and navigating a device menu.**

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