PRODUCT INSERT

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NOTE: *The following information is in addition to the information contained in the Power Lift Recliner Owner's Manual provided with the power lift recliner. It explains how the Patent-Pending FirstUp Technology affects the speed at which your power lift recliner operates. Please retain this information with the Power Lift Recliner Owner's Manual and contact your authorized Pride Provider with any questions.*

If your power lift recliner was purchased with factory-installed Patent-Pending FirstUp Technology, this feature reduces the time required to raise your power lift recliner from a seated or reclined position between 15 to 19 seconds, depending on user weight. Patent-Pending FirstUp Technology does not affect the time required to lower or recline your power lift recliner.

NOTE: *If the power lift recliner weight capacity is exceeded, the power lift recliner may operate at normal speed.*

NOTE: *Using the battery backup will cause the power lift recliner to operate at normal speed.*

NOTE: *Patent-Pending FirstUp Technology has been designed specifically to operate with specialized manufacturer-provided components.*

NOTE: *When reordering components, use the following part numbers:*

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