

Backpacker PLUS

Base weight lift	155 lbs. ¹
Scooter/power chair overall weight	325 lbs.
Maximum product width	27"
Maximum travel distance from the ground	26"
Motor	Sealed 12 volt DC
Overall platform length	44" (41" usable surface)
Maximum product wheelbase	35"
Warranty	3-year transferable
Optional battery pack	12 volt, on-board
Optional charger	110 volt

¹ Weight does not include any accessories
² Vehicle dependent

Can easily lift the following current Pride® Power Chairs²

- Jazzy® 600
- Jazzy® 614/614HD
- Jazzy® 1103 Ultra
- Jazzy® 1113 ATS
- Jazzy Select®
- Jazzy Select® 14
- Jazzy Select® 6
- Jazzy Select® 6 Ultra
- Jazzy Select® Traveller
- Jet™ 3 Ultra

Can easily lift the following current Quantum® Power Chairs²

- Quantum® 600/600 XL
- Quantum® 610
- Quantum® 614
- Quantum® 1103 Ultra
- Quantum® 1121
- Quantum® R4000

Can easily lift the following current Pride® Scooters²

- Celebrity® X 3 & 4
- Go-Go® Ultra X 3 & 4
- Go-Go Elite Traveller®/ LR 3 & 4
- Go-Go Elite Traveller® Plus 3 & 4
- Go-Go Elite Traveller® Plus HD 3 & 4
- Legend® 3 & 4
- Maxima 3 & 4
- Rally® 3 & 4
- Revo® 3 & 4
- Sundancer®
- Victory® 10 3 & 4

¹ Weight does not include any accessories
² Vehicle dependent

Backpacker PLUS



Battery Pack

NEW, OPTIONAL MULTI-PURPOSE BATTERY BOX



- Industry first
- Industry exclusive
- Charges the mobility device from the lift as you drive



Backpacker PLUS

Easy to install and use, the Silver Star® Backpacker® Plus Interior Platform Lift features a pre-assembled single piece platform that you can drive your scooter or power chair on or off from either direction. Enhanced construction and a fresh new look give the Backpacker Plus the quality you expect and deserve from Silver Star Lifts.

- The quality construction and design you expect from Pride® Mobility Products
- 325 lbs. lifting capacity
- Easy-to-use 2-button hand control
- 3-year transferable warranty
- Can lift a variety of Pride® Scooters and Power Chairs and other manufacturers' mobility devices

Note: The Backpacker® Plus can lift all of the listed power chairs with contoured, Synergy® or Specialty seating. Please call Customer Service when requesting information about lifting with a power positioning system. Can fit in most mini-vans with little or no vehicle modification, full-size vans require modification.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



All black battery box, bracket, and new matte finish



One piece, fully assembled drive-on/drive-off platform



Single piece horizontal drive rail: Minimum welding points for noise reduction, increased rigidity, and more fluid movement



Removable, front drive panel for easier serviceability



4" more usable space



New, angled, low-profile platform



New, adjustable wheel chocks on the platform surface