

# GO-GO

ELITE TRAVELLER® PLUS



*Go-Go Elite Traveller® Plus 3-wheel*

---

The **Go-Go Elite Traveller® Plus** was designed to bring travel scooters to a greater range of people than ever before. An increase in length and width and a weight capacity of 325 lbs. combine with a wraparound delta tiller to allow both the larger individual and those with limited dexterity to experience the benefits of a travel specific mobility product for the first time. One hand feather-touch disassembly makes the Go-Go Elite Traveller Plus the easiest travel scooter to take along anywhere.

---



*Go-Go Elite Traveller® Plus 4-wheel*

Available as a 3 & 4-wheel scooter



**Feather-touch disassembly.**

**Pride**<sup>®</sup>  
Mobility Products Co.  
— LIVE YOUR BEST —



Max. weight capacity	325 lbs.
Maximum speed <sup>1</sup>	Up to 4.5 mph
Ground clearance <sup>2</sup>	3-wheel: 3" 4-wheel: 3"
Turning radius <sup>2</sup>	3-wheel: 34.5" 4-wheel: 45.75"
Overall length <sup>2</sup>	3-wheel: 40" 4-wheel: 40.25"
Overall width <sup>2</sup>	3-wheel: 21.125" 4-wheel: 21.125"
<b>Tires</b>	
Front	3-wheel: 2" x 8" solid 4-wheel: 2" x 8" solid
Rear	3-wheel: 3" x 9" solid 4-wheel: 3" x 9" solid
Anti-tips	Rear
Standard seating	Compact, foldable, black vinyl
Suspension	NA
Battery charger	Off-board, 2 amp
Per-charge range <sup>1,5</sup>	Up to 12 miles
Battery requirements <sup>3</sup>	(2) 12V deep-cycle
Compatible batteries	17AH
<b>Weight</b>	
Battery weight <sup>4</sup>	17AH battery pack: 29 lbs.
Weight w/o battery pack	3-wheel: 84 lbs. 4-wheel: 87 lbs.
Heaviest piece	3-wheel: 32.5 lbs. (rear section) 4-wheel: 35 lbs. (front section)

**Colors**



Red



Blue

2 sets of exciting colors with every scooter!


**Features**

- Feather-touch disassembly with one hand
- Delta tiller with wraparound handles
- Pride's exclusive black, non-scuffing tires
- Standard 17 AH battery pack
- 9" rear wheels
- 325 lbs. weight capacity
- Auto-connecting front to rear cable
- Includes sets of easily changeable colored panels - Red and Blue
- Front frame-mounted seat post offers maximum stability
- Frame design easily disassembles into 5 super lightweight pieces for convenient transport and storage
- Modular design for easy serviceability
- Microprocessor-based controller offers optimal power management and added safety features
- Front tiller-mounted basket
- Removable footboard basket for additional storage (4-wheel version only)
- Convenient off-board dual voltage charger can charge battery pack on-board or off-board

**Pride™ Lifts & Ramps**

These Pride™ Lifts are compatible with the Go-Go Elite Traveller® Plus.

- Backpacker® Plus
- Outlander
- SilverBoom 250
- Power Tote Mini (3-wheel only)



SilverBoom 250



Backpacker® MV

St. Catharines, Ontario • 888-570-1113 • www.pridemobility.com

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.  
<sup>2</sup> Due to manufacturing tolerances and continued product improvement, this specification is subject to a variance of +/- 3%.  
<sup>3</sup> AGM or gel cell type recommended.  
<sup>4</sup> Battery weight may vary based on manufacturer.  
<sup>5</sup> Tested in accordance with ANSI/RENSA, WC Vol 2, Section 4, and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.

