

# Jazzy® 1113 ATS

The Jazzy® 1113 ATS offers everything you need—all of the “big chair” options and accessories, even rehab capabilities—on the most compact frame available. You’ll take on tight spaces and narrow hallways without sacrificing a strong outdoor performance never before seen in a chair of this size. For sharp ramp transitions as well as a smooth and stable ride even over exceptionally rough and uneven terrain, the Jazzy 1113 ATS is equipped with Active-Trac® Suspension.



Jazzy® 1113 ATS  
shown in Candy Apple Red

St. Catharines, Ontario ■ 888-570-1113  
Visit us on the web: [www.pridemobility.com](http://www.pridemobility.com)

**Pride**®  
Mobility Products Co.

## SPECIFICATIONS

<b>Weight capacity</b>	300 lbs.
<b>Maximum speed</b>	Up to 4 mph
<b>Ground clearance</b>	2.5"
<b>Turning radius</b>	20.5"
<b>Overall length</b>	36"
<b>Overall width</b>	23"
<b>Drive Wheels</b>	10"
<b>Front anti-tips</b>	6" solid
<b>Rear casters</b>	8" solid
<b>Suspension</b>	Active-Trac® Suspension
<b>Drive train</b>	Two-motor, Mid-Wheel Drive
<b>Braking system</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard electronics</b>	50 amp, PG VSI controller
<b>Specialty controls</b>	Yes
<b>Battery charger</b>	On-board, 4 amp (standard)
<b>Per-charge range</b>	Up to 15 miles
<b>Battery requirements</b>	(2) 12 volt, deep cycle
<b>Battery size</b>	U-1
<b>Battery weight</b>	24.5 lbs. each
<b>Base weight</b>	77 lbs.
<b>Standard seat weight</b>	37 lbs. (medium-back)
<b>Max. seat dimensions</b>	20" x 20"
<b>Warranty</b>	Lifetime limited warranty on frame; 2-year limited warranty on electronics; 18-month limited warranty on drive motors

## COLOR CHOICES



Candy Apple  
Red



Viper Blue



Onyx Black



Champagne

## FEATURES

- Powerful two-motor, mid-wheel drive design offers outstanding maneuverability
- Low center of gravity design delivers a very smooth and stable ride
- 10" knobby drive wheels and optimum weight distribution provide superior traction
- On-board battery charger
- Fully programmable controller

## OPTIONS


- TRU-Balance® Power Tilt <sup>1</sup>
- TRU-Balance® Power Recline <sup>1</sup>
- Synergy® seating system
- Swing-away joystick mount
- Solid seat pan
- Cane/crutch holder
- Cup holder
- Walker holder
- Weather cover
- Safety flag
- Rear basket
- Off-board charger
- Oxygen tank holder
- Solid seat pan

<sup>1</sup>) 200 lbs. weight capacity with TRU-Balance® Systems

## RECOMMENDED LIFTS

**Pride SilverStar**  
**Lifts & Ramps** 

- Backpacker® Plus
- Backpacker® MV
- Full Platform
- Silver Boom 300



**Pride**®  
Mobility Products Corp.

St. Catharines, Ontario • 888-570-1113  
Visit us on the web: [www.pridemobility.com](http://www.pridemobility.com)

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

\* Special order seats or non-standard dimensions will not be fold-flat and depth adjustable.