

Jazzy®

1170 XL PLUS

The Jazzy® 1170XL Plus is a stylish, high performance power chair that combines excellent indoor maneuverability with a top speed up to 5 mph and a weight capacity of up to 500 lbs. It is equipped standard with Active-Trac® Suspension which provides for great outdoor performance and rough terrain handling. The Jazzy 1170XL Plus is the ultimate power chair for the active outdoor user who still requires tight-quarter indoor maneuverability.



Jazzy® 1170 XL Plus
shown in *Candy Apple Red*

St. Catharines, Ontario ■ 888-570-1113
Visit us on the web: www.pridemobility.com

Pride®
Mobility Products Co.

SPECIFICATIONS

Weight capacity	500 lbs.
Maximum speed	Up to 5 mph
Ground clearance	4"
Turning radius	26.75"
Overall length	45.25"
Overall width	30"
Drive Wheels	16" solid (required)
Front anti-tips	8" solid
Rear casters	9" solid or pneumatic
Suspension	Active-Trac® Suspension
Drive train	Two-motor, Mid-Wheel Drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	70 amp, PG VSI controller
Specialty controls	Yes
Battery charger	Off-board, 8 amp (standard)
Per-charge range	Up to 15 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	Group 24
Battery weight	53.5 lbs. each
Base weight	149.5 lbs.
Standard seat weight	42.5 lbs.
Warranty	Lifetime limited warranty on frame; 2-year limited warranty on electronics; 18-month limited warranty on drive motors

COLOR CHOICES



Candy Apple
Red



Viper Blue



Onyx Black



Champagne

FEATURES

- 500 lbs. weight capacity
- Powerful, patented, two-motor mid-wheel drive design provides superb power and acceleration
- Active-Trac® and rear caster suspension for outstanding outdoor performance
- Adjustable front anti-tips to meet your preferences and terrain challenges
- On-board battery charger
- Battery-conserving sleep mode available with Remote Plus controller

OPTIONS

- Contoured high-back seat
- Cane/crutch holder
- Lap belt
- Heavy-duty dust cover
- Cup holder
- Contoured medium/high-back with solid seat pan
- Rear basket
- Oxygen tank holder
- Walker holder
- Solid tires
- Safety flag
- Swing-away joystick mount

RECOMMENDED LIFTS

Silver Star®
Lifts & Ramps

- Backpacker® Plus
- Backpacker® MV
- Full Platform
- Silver Boom 300



PRIDE
PERFORMANCE MOBILITY
— LIVE YOUR BEST —

St. Catharines, Ontario • 888-570-1113
Visit us on the web: www.pridemobility.com

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

* Special order seats or non-standard dimensions will not be fold-flat and depth adjustable.