

Jazzy® 600

The Jazzy® 600 provides excellent performance and maximum stability. Its compact design and standard ATX Suspension allows it to excel indoors and out making it a superb all-around power chair choice.



Jazzy® 600
shown in *Candy Apple Red*

St. Catharines, Ontario ■ 888-570-1113
Visit us on the web: www.pridemobility.com

Pride®
Mobility Products Co.

SPECIFICATIONS

Weight capacity	300 lbs.
Maximum speed	Up to 5 mph
Ground clearance	3"
Turning radius	20"
Overall length	36.5" without front riggings
Overall width	23"
Wheels	Drive: 14"
Front casters	5" solid OMNI-casters
Rear casters	6" solid OMNI-casters
Suspension	Active-Trac® ATX Suspension
Drive train	Two-motor, Mid-Wheel 6®
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	60 amp, PG VR2
Specialty controls	N/A
Battery charger	Off-board, 5 amp
Per-charge range	Up to 20 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	NF-22
Battery weight	38 lbs. each
Base weight	127 lbs.
Standard seat weight	37 lbs. (medium-back)
Max. seat dimensions	20" x 20"
Seat-to-floor range (high-back seat)	21" - 24"
Warranty	Lifetime limited warranty on frame; 2-year limited warranty on electronics; 18-month limited warranty on drive motors

FEATURES

- OMNI-Casters (nylon, spherical-shaped casters) on front and rear to prevent wheel hang-ups
- Side-mounted, easily accessible freewheel levers
- Active-Trac® ATX Suspension incorporates front OMNI-Casters for enhanced performance over more varied terrain
- Easy front access to batteries

OPTIONS

- Solid seat pan
- Elevating leg rests
- Solid tires
- Oxygen tank holder
- Cup holder
- Cane/crutch holder
- Swing-away footrests
- Walker holder
- Swing-away joystick
- Contoured high-back seat

COMPATIBLE LIFTS



- Backpacker® Plus
- Backpacker® MV
- Full Platform
- Silver Boom 300

COLOR CHOICES



Candy Apple Red Viper Blue Onyx Black Champagne