

# Jazzy®

614 | 614HD

The Jazzy® 614 and 614HD deliver outstanding performance and superb maneuverability in a powerful, stylish package. Enhanced in-line motors, advanced ATX suspension and upgraded electronics make both the Jazzy 614 and 614HD innovative and exceptional overall values.



Jazzy® 614  
shown in Jazzy Red


St. Catharines, Ontario • 888-570-1113  
Visit us on the web: [www.pridemobility.com](http://www.pridemobility.com)

**Pride®**  
Mobility Products Co.

## SPECIFICATIONS

|                             |   |
|-----------------------------|---|
| <b>Weight capacity</b>      | 300 lbs. (614 HD: 450 lbs.)   |
| <b>Maximum speed</b>        | Up to 4.25 mph  |
| <b>Ground clearance</b>     | 2.25"   |
| <b>Turning radius</b>       | 20"   |
| <b>Overall length</b>       | 35.5" without front riggings  |
| <b>Overall width</b>        | 25.125"   |
| <b>Wheels</b>               | Drive: 14"  |
| <b>Front casters</b>        | 6" solid  |
| <b>Rear casters</b>         | 6" solid  |
| <b>Suspension</b>           | Active-Trac® ATX Suspension   |
| <b>Drive train</b>          | Two-motor, in-line, Mid-Wheel 6®  |
| <b>Braking system</b>       | Intelligent braking<br>(electronic, regenerative disc brakes)   |
| <b>Standard electronics</b> | 50 amp, Dynamic controller  |
| <b>Battery charger</b>      | Off-board, 5 amp (standard)   |
| <b>Per-charge range</b>     | Up to 16 miles (614 HD: up to 11 miles)   |
| <b>Battery requirements</b> | (2) 12 volt, deep cycle   |
| <b>Battery size</b>         | NF-22   |
| <b>Battery weight</b>       | 38 lbs. each  |
| <b>Base weight</b>          | 121.5 lbs.  |
| <b>Standard seat weight</b> | 37 lbs. (medium-back)<br>42.5 lbs. (high-back on 614HD)   |
| <b>Max. seat dimensions</b> | 22" x 22" (614 HD: 24" x 24")   |
| <b>Warranty</b>             | 5-year limited warranty on frame;<br>1-year limited warranty on electronics;<br>1-year limited warranty on drive motors |

## COLOR CHOICES



Jazzy Red



Jazzy Blue

## FEATURES

- New, dual enhanced in-line motors for greater efficiency, torque, range and performance
- 450 lbs. weight capacity with HD model
- 6" casters on front & rear for maximum stability
- Active-Trac® ATX Suspension for enhanced performance over more varied terrain
- Easy rear access to batteries
- Seating up to 24" W x 24" D with HD model
- Tool-less shroud removal
- Easy serviceability
- Depth adjustable foot platform
- Fold-flat, depth adjustable seat \*  
(Standard on J614HD; Optional on J614)

## OPTIONS

- Solid seat pan
- Cup holder
- Swing-away footrests
- Elevating leg rests
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

## RECOMMENDED LIFTS

**Pride SilverStar®**  
**Lifts & Ramps** 

- Backpacker® MV
- Backpacker® Plus
- Full platform
- Silver Boom 300



**Pride**<sup>®</sup>  
Mobility Products Co.

St. Catharines, Ontario • 888-570-1113  
Visit us on the web: [www.pridemobility.com](http://www.pridemobility.com)

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

\* Special order seats or non-standard dimensions will not be fold-flat and depth adjustable.