

# Jazzy<sup>®</sup>

## SELECT 14

The Jazzy Select<sup>®</sup> 14 Series delivers the advantages of patented mid-wheel drive for maneuvering in tight, compact spaces. The Jazzy Select 14 is loaded with performance and convenience features like a powerful 60 amp controller, larger batteries for greater per-charge range, and an on-board battery charger.



Jazzy Select<sup>®</sup> 14  
shown in Jazzy Red

St. Catharines, Ontario • 888-570-1113  
Visit us on the web: [www.pridemobility.com](http://www.pridemobility.com)

**Pride<sup>®</sup>**  
Mobility Products Co.

## SPECIFICATIONS

<b>Weight capacity</b>	300 lbs.
<b>Maximum speed</b>	Up to 5 mph
<b>Ground clearance</b>	3"
<b>Turning radius</b>	21"
<b>Overall length</b>	38" without front riggings
<b>Overall width</b>	24.75"
<b>Drive wheels</b>	14"
<b>Front casters</b>	6" solid
<b>Rear casters</b>	8" solid
<b>Suspension</b>	Limited
<b>Drivetrain</b>	Two-motor, mid-wheel drive
<b>Braking system</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard electronics</b>	60 amp controller
<b>Specialty controls</b>	N/A
<b>Battery charger</b>	Off-board, 5 amp
<b>Per-charge range</b>	Up to 15 miles
<b>Battery requirements</b>	(2) 12 volt, deep cycle
<b>Battery size</b>	NF-22
<b>Battery weight</b>	38 lbs. each
<b>Base weight</b>	105 lbs.
<b>Seat weight</b>	37 lbs. (medium-back)
<b>Max. seat dimensions</b>	20" x 20"
<b>Warranty</b>	5-year limited warranty on frame; 1-year limited warranty on electronics; 1-year limited warranty on drive motors

## COLOR CHOICES



Jazzy Red



Jazzy Blue

St. Catharines, Ontario • 888-570-1113  
Visit us on the web: [www.pridemobility.com](http://www.pridemobility.com)

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

## FEATURES

- Patented mid-wheel drive design for tight space maneuverability
- 14" knobby tires and optimum weight distribution provide superior traction
- NF-22 batteries for maximum range
- 5 amp, off-board charger

## OPTIONS

- Depth-adjustable high-back seat with headrest
- High-back seat with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

## RECOMMENDED LIFTS

**Silver Star®**  
**Lifts & Ramps**

- Backpacker® AVP
- Backpacker® Plus
- Backpacker® MV
- Full Platform
- Silver Boom 300



**PRIDE**  
PERFORMANCE MOBILITY

— LIVE YOUR BEST® —