

Jazzy[®] SELECT

The easily maneuverable, restyled Jazzy Select[®] is equipped with performance features like Active-Trac[®] Suspension and in-line motor technology for enhanced efficiency, torque and range for high-performance. With features like a Contour high-back seat with a headrest, full length armrests and black, non-scuffing tires with silver rims, the Jazzy Select is a better overall value than ever before.



Jazzy Select[®] shown in Jazzy Blue with gray vinyl medium-back seat, PG GC 3 controller and foot platform

NEW
FOR 2010

St. Catharines, Ontario • 888-570-1113
Visit us on the web: www.pridemobility.com

Pride[®]
Mobility Products Co.

SPECIFICATIONS

Weight capacity	300 lbs.
Maximum speed	Up to 4 mph
Ground clearance	1.875" at center of frame
Turning radius	20"
Overall length	35.5" without front riggings
Overall width	23.75"
Drive wheels	10"
Front anti-tips	5" solid
Rear casters	6" solid
Suspension	Active-Trac®
Drivetrain	Two-motor, in-line mid-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	40 amp, PG GC 3 controller
Specialty controls	N/A
Battery charger	Off-board, 3 amp
Per-charge range	Up to 15 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	U-1
Battery weight	24.5 lbs. each
Base weight	99 lbs.
Standard seat weight	37 lbs. (medium-back, contour seat)
Max. seat dimensions	20" x 20"
Warranty	5-year limited warranty on frame; 1-year limited warranty on electronics; 1-year limited warranty on drive motors

COLOR CHOICES



Jazzy Red



Jazzy Blue

St. Catharines, Ontario • 888-570-1113
Visit us on the web: www.pridemobility.com

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

FEATURES

- New black, non-scuffing tires with silver rims
- Larger foot platform
- 40 amp, PG GC 3 controller
- Built with ease of service in mind
- Easy side access to batteries
- Toolless front panel access to electronics to minimize service time
- Active-Trac® Suspension
- Fold-flat, depth-adjustable high-back
- Medium-back seat

OPTIONS

- High-back seat
- Medium-back seat with solid seat pan
- Limited-recline high-back seat
- Limited-recline high-back seat with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

RECOMMENDED LIFTS

Silver Star®
Lifts & Ramps

- Backpacker® Plus
- Full Platform
- Backpacker® MV
- Silver Boom 300



PRIDE
PERFORMANCE MOBILITY

— LIVE YOUR BEST® —