

JET[®] 3


ULTRA

The Jet[®] 3 Ultra delivers a broad range of standard performance features like Active-Trac[®] Suspension and a new, compact 40 amp, PG GC 2 controller on a very maneuverable frame. Plus, the style is fresh, exciting and contemporary. These things combine with a truly economical price to make the Jet 3 Ultra an excellent overall value in personal mobility.



	Jet® 3 Ultra
Weight capacity:	300 lbs.
Maximum speed:	Up to 4 mph
Ground clearance:	2.75" at frame
Turning radius:	20"
Overall length:	35.5" without front rigging
Overall width:	22.5"
Wheels:	Drive: 10" solid
Front Anti-tips:	6" solid
Rear Casters:	8" solid
Suspension:	Active-Trac®
Drive train:	Two-motor, mid-wheel drive
Braking system:	Intelligent braking (electronic, regenerative disc brakes)
Std. Electronics:	40 amp, PG GC 2 controller
Specialty controls:	N/A
Battery charger:	Off-board, 3 amp
Per-charge range:	Up to 11 miles
Battery requirements:	(2) 12 volt, deep cycle
Battery size:	U-1
Transport weight	
Battery weight:	24.5 lbs. each
Base weight:	92 lbs.
Std. seat weight:	37 lbs. (medium-back)
Std. Seat dimensions:	18" x 18"
Warranty:	5-year limited warranty on frame; 1-year limited warranty on electronics; 1-year limited warranty on drive motors

Colors:



Jet Red



Jet Blue



FEATURES

- Compact size for superior handling in tight quarters
- Standard Active-Trac® Suspension
- 300 lbs. weight capacity
- Utilizes trapeze bar seat mount for quick adjustability and easy interchange of a broad range of seats
- Convenient front battery access
- Curved rear caster beam provides excellent stability
- Maximum speed up to 4 mph
- Up to 15 miles of range per charge
- 20" turning radius

OPTIONS

- High-back seat
- High-back seat with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Swing-away footrests
- Elevating leg rests



COMPATIBLE LIFTS

- Backpacker® Plus
- Backpacker® MV
- Backpacker® AVP
- Full Platform
- Silver Boom



St. Catharines, Ontario • 888-570-1113
Visit us on the web: www.pridemobility.com

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %.