

GOCHAIR®

Pride®
Mobility Products Corp.
— LIVE YOUR BEST™ —

The first travel chair that combines super portability, maneuverability and style at a low, low cost!



Specifications

Drive Wheels	8" (solid)
Front Casters	5" (solid)
Rear Anti-tips	3" (solid)
Maximum Speed	3.5 mph
Ground Clearance	1.75"
Turning Radius	26.75" with foot platform
Length	29" without foot platform
Width	19"
Seat Type	Foldable, compact, lightweight
Batteries	(2) 12V, 12 amp
Range	Up to 8 miles
Charger	Off-board, 2 amp
Electronics	40 amp, PG GC 2 (non-programmable)
Weight Capacity	250 lbs.
Base Weight	57 lbs. (front & rear)
Standard Seat Weight	25 lbs.
Battery Pack Weight	20 lbs.
Suspension	Limited
Heaviest Piece (rear)	34.5 lbs.

Portability

- Easy-to-remove battery pack
- Conveniently disassembles into four manageable pieces
- Fits in most vehicle trunks

Performance

- Top speed of 3.5 mph and enough power to get you where you want to go
- Twin motors for traction and precision handling

Maneuverability

- Compact design for a tight turning radius

Colors



Red



Blue

(US) 800-800-8586 • (Canada) 888-570-1113 • www.pridemobility.com

Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. This product is not coded by DME PDAC or does not meet criteria: K0899.

PRIDE®
PERFORMANCE MOBILITY
— LIVE YOUR BEST™ —