

# GO-GO

## ELITE TRAVELLER® PLUS AND ELITE TRAVELLER® PLUS HD



*Go-Go Elite Traveller® Plus, 3-Wheel*

---

The **Go-Go Elite Traveller® Plus** was designed to bring travel scooters to a greater range of people than ever before. An increase in length and width, and a weight capacity of 300 lbs. (325 lbs. w/HD version) combine with wraparound delta tiller to allow both the larger individual and those with limited dexterity to experience the benefits of a travel-specific mobility product for the first time. One hand feather-touch disassembly makes the Go-Go Elite Traveller Plus the easiest travel scooter to take along anywhere.

---



*Go-Go Elite Traveller® Plus,  
4-Wheel*

**Available as a 3 & 4-Wheel scooter**



**Feather-touch disassembly.**

**Pride®**  
Mobility Products Corp.  
— LIVE YOUR BEST™ —



ELITE TRAVELLER PLUS  
AND

ELITE TRAVELLER PLUS HD

Max. weight capacity	300 lbs./325 lbs. w/HD version
Maximum speed <sup>1</sup>	Up to 4.5 mph
Ground clearance <sup>2</sup>	3-Wheel: 3" 4-Wheel: 3"
Turning radius <sup>2</sup>	3-Wheel: 34.5" 4-Wheel: 45.75"
Overall length <sup>2</sup>	3-Wheel: 40" 4-Wheel: 40.25"
Overall width <sup>2</sup>	3-Wheel: 21.125" 4-Wheel: 21.125"
<b>Tires</b>	
Front	3-Wheel: 2" x 8" solid 4-Wheel: 2" x 8" solid
Rear	3-Wheel: 3" x 9" solid 4-Wheel: 3" x 9" solid
Anti-tips	Rear
Standard seating	Compact, foldable, black vinyl
Suspension	NA
Battery charger	Off-board, 2 amp
Per-charge range <sup>1,5</sup>	Up to 12 miles
Battery requirements <sup>3</sup>	(2) 12V deep-cycle
Compatible batteries	18AH
<b>Weight</b>	
Battery weight <sup>4</sup>	18AH battery pack: 29 lbs.
Weight w/o battery pack	3-Wheel: 84 lbs. 4-Wheel: 87 lbs.
Heaviest piece	3-Wheel: 32.5 lbs. (rear section) 4-Wheel: 35 lbs. (front section)

**Standard Colors**



Red



Blue

2 sets of exciting colors with every scooter!


**Features**

- Feather-touch disassembly with one hand
- Delta tiller with wraparound handles
- Pride's exclusive black non-scuffing tires
- Standard 18 AH battery pack
- 9" rear wheels
- 325 lbs. weight capacity w/HD version
- Auto-connecting front to rear cable
- Includes sets of easily changeable colored panels - Red and Blue
- Front frame-mounted seat post offers maximum stability
- All new frame design easily disassembles into 5 super lightweight pieces for convenient transport and storage
- Modular design for easy serviceability
- Microprocessor-based controller offers optimal power management and added safety features
- Front tiller-mounted basket
- Removable footboard basket for additional storage (4-Wheel version only)
- Convenient off-board dual voltage charger can charge battery pack on-board or off-board

**Silver Star Lifts & Ramps**

These Silver Star® Lifts are compatible with the Go-Go Elite Traveller® Plus.

- Backpacker® Plus
- Full Platform



Silver Boom 300



Backpacker® MV

<sup>1</sup>Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.

<sup>2</sup>Due to manufacturing tolerances and continued product improvement, this specification is subject to a variance of +/- 3%.

<sup>3</sup>AGM or gel cell type recommended.

<sup>4</sup>Battery weight may vary based on manufacturer.

<sup>5</sup>Tested in accordance with ANSI/RENSA, WC Vol 2, Section 4, and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.

(US) 800-800-8586 • (Canada) 888-570-1113 • [www.pridemobility.com](http://www.pridemobility.com)