

FLASH


REAR-WHEEL DRIVE POWER CHAIR



The Flash's compact rear-wheel drive design and shorter wheelbase provide tight-space maneuverability with a turning radius of only 22.25 inches! Additionally, the Flash offers many exciting standard features and sporty styling all in a very economical package.



182 Susquehanna Avenue
Exeter, PA 18643
(U.S.) 1-800-800-8586
(Canada) 1-888-570-1113
www.pridemobility.com



SPECIFICATIONS

Weight capacity:	300 lbs.
Maximum speed:	Up to 5 mph
Ground clearance:	3"
Turning radius:	22.25"
Overall length:	29"
Overall width:	24"
Wheels:	Drive: 10" solid or pneumatic Front casters: 6" solid Anti-tips: 3" solid
Suspension:	Limited
Drive train:	Two-motor, rear-wheel drive
Braking system:	Intelligent braking (electronic, regenerative disc brakes)
Std. Electronics:	50 amp, PG VSI controller
Specialty controls:	No
Battery charger:	On-board, 4 amp (standard)
Per-charge range:	Up to 20 miles
Battery requirements:	(2) 12 volt, deep cycle
Battery size:	32 AH, U-1
Transport weight	
Battery weight:	24.5 lbs. each
Base weight:	89.5 lbs.
Seating:	Gray vinyl contoured medium-back
Seat weight:	33 lbs.
Seat dimensions:	18" W x 18" D (max. medium-back dim.)
Warranty:	Lifetime limited warranty on frame 1-year limited warranty on electronics 1-year limited warranty on drive motors



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

FEATURES

- Tight-space maneuverability with only a 22.25" turning radius
- Up to 300 lbs. weight capacity
- Swing-away joystick
- 180° swivel seat
- On-board battery charger with easy access to rear charger panel for A/C receptacle and ammeter
- Articulating front caster beam for maximum stability
- Extended freewheel release handles for easy operation
- Seat-to-floor height as low as 20.25"



Pride
Mobility Products Corp.