

Jazzy® 1103 ULTRA

The Jazzy® 1103 Ultra brings sporty style and high performance features, like a standard power elevating seat and Active-Trac® Suspension, together into one perfectly balanced power chair. Yet, even fully loaded, the Jazzy 1103 Ultra remains highly compact and maneuverable.



Jazzy® 1103 Ultra
shown in *Candy Apple Red*

(US) 800-800-8586 ▪ (Canada) 888-570-1113
Visit us on the web: www.pridemobility.com

Pride®
Mobility Products Corp.

SPECIFICATIONS

Weight capacity	300 lbs.
Maximum speed	Up to 4 mph
Ground clearance	2.5"
Turning radius	20.75"
Overall length	35.5"
Overall width	23"
Drive Wheels	10"
Front anti-tips	5" solid
Rear casters	8" solid
Suspension	Active-Trac® Suspension
Drive train	Two-motor, Mid-Wheel Drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	50 amp, 6-key PG VSI controller
Specialty controls	No
Battery charger	On-board, 4 amp (standard)
Per-charge range	Up to 15 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	U-1
Battery weight	24.5 lbs. each
Base weight	104 lbs.
Standard seat weight	37 lbs. (medium-back)
Max. seat dimensions	20" x 20"
Warranty	Lifetime limited warranty on frame; 2-year limited warranty on electronics; 18-month limited warranty on drive motors


COLOR CHOICES



Candy Apple
Red



Viper Blue



Onyx Black



Champagne

MODELS

Jazzy 1103 Ultra-C
Jazzy 1103 Ultra-SS

FEATURES

- 300 lbs. weight capacity
- Standard power elevating seat provides greater access to objects out of reach
- Powerful two-motor, patented mid-wheel drive design offers outstanding maneuverability
- Low center of gravity design delivers a very smooth and stable ride
- 10" knobby drive wheels and optimum weight distribution provide superior traction
- Speed range adjustable up to 4 mph
- Standard lap belt
- On-board battery charger


OPTIONS

- Contoured high-back seat
- Swing-away joystick mount
- Solid seat pan
- Cane/crutch holder
- Cup holder
- Walker holder
- Weather cover
- Safety flag
- Rear basket
- Off-board charger
- Oxygen tank holder
- Solid seat pan

RECOMMENDED LIFTS

Pride SilverStar®
Lifts & Ramps

- Backpacker® Plus
- Backpacker® MV
- Full platform
- Silver Boom 300



Pride®
Mobility Products Corp.

(US) 800-800-8586 ▪ (Canada) 888-570-1113
Visit us on the web: www.pridemobility.com

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

* Special order seats or non-standard dimensions will not be fold-flat and depth adjustable.