

Jazzy®

1133

The sporty Jazzy® 1133 delivers high-performance, hot style, rehab capability, and patented features, like Active-Trac® Suspension, in a full-size, yet easily transportable package that conveniently disassembles into five manageable base pieces. All these features combine with true rigid-frame stability to bring innovation to the world of full-sized transportable power chairs. Simply put, the Jazzy 1133 is the perfect full-size, take-along power chair choice for people on the go.



DISASSEMBLY

Quickly and easily disassembles



VERSATILITY

Rehab capabilities



PERFORMANCE

Active-Trac® Suspension for excellent outdoor performance



STABILITY

Full-sized rigid frame stability

Pride
Mobility Products Corp.


800-800-8586 (US)
888-570-1113 (Canada)
pridemobility.com

Weight capacity:	300 lbs.
Maximum speed:	Up to 5 mph
Ground clearance:	3.5"
Turning radius:	21.5"
Overall length:	41"
Overall width:	25.5"
Wheels:	Drive: 14" pneumatic (standard)
	Casters: 8" solid
	Anti-tips: 6" solid
Suspension:	Active-Trac®
Drive train:	Two-motor, mid-wheel drive
Braking system:	Intelligent braking (electronic, regenerative disc brakes)
Std. Electronics:	50 amp, PG VSI controller
Specialty controls:	N/A
Battery charger:	On-board, 4 amp (standard)
Per-charge range:	Up to 25 miles
Battery requirements:	(2) 12 volt deep cycle
Battery size:	NF-22
Transport weight	
Battery weight:	38 lbs. each
Base weight:	95 lbs.
Std. seat weight:	37 lbs. (medium-back)
Seat dimensions:	20" x 22" (max. contoured seating)
Warranty:	Lifetime limited warranty on frame 2-year limited warranty on electronics 18-months limited warranty on drive motors

Colors:



Easy disassembly for quick and convenient transport



FEATURES

- Easily disassembles for transport
- Powerful two-motor, mid-wheel drive design offers outstanding maneuverability
- Active-Trac® Suspension for outstanding performance
- Adjustable front anti-tips for personal preferences and terrain challenges
- Speed range adjustable up to 5 mph
- On-board battery charger
- Fully programmable controller
- Easy rear battery access

OPTIONS

- Synergy® seating
- Contoured high-back seat
- Lap belt
- Cup holder
- Dust cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Flat-free inserts
- Adjustable-height armrests
- Flat-free tires
- Saddlebag

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.