

Jazzy®

1143 ULTRA

ENHANCED OUTDOOR PERFORMANCE IN A COMPACT DESIGN

The Jazzy 1143 Ultra combines sporty style and rugged outdoor performance with a standard power elevating seat. The 12" drive wheels offer enhanced outdoor capability, while its mid-wheel drive maneuverability and smaller frame equal a tight turning radius for easy indoor use. All of these features add up to a power chair that's ready for your active lifestyle.



Jazzy® 1143 Ultra shown in Candy Apple Red with power elevating gray contoured high-back seat, foot platform and VSI controller.

Pride
Mobility Products Corp.

800-800-8586 (US)
888-570-1113 (Canada)
pridemobility.com



PERFORMANCE

12" drive wheels and Active-Trac® Suspension for superb outdoor performance



INNOVATION

Standard power elevating seat provides greater access to your environment



STYLE

Sleek sporty look



MANEUVERABILITY

Compact design and tight turning radius are perfect for small spaces

SPECIFICATIONS

Weight capacity:	300 lbs.
Maximum speed:	Up to 4.5 mph
Ground clearance:	3"
Turning radius:	21.25"
Overall length:	36"
Overall width:	24.75"
Wheels:	Drive: 12" pneumatic Casters: 8" solid Anti-tips: 6" solid
Suspension:	Active-Trac®
Drive train:	Two-motor, mid-wheel drive
Braking system:	Intelligent braking (electronic, regenerative disc brakes)
Std. Electronics:	50 amp, 6-key VSI
Specialty controls:	N/A
Battery charger:	On-board, 4 amp
Per-charge range:	Up to 25 miles
Battery requirements:	(2) 12 volt deep cycle
Battery size:	U-1
Transport weight:	
Battery weight:	24.5 lbs. each
Base weight:	112.5 lbs.
Std. seat weight:	37 lbs. (medium-back)
Seat dimensions:	20" x 20" (max. contoured seating)
Warranty:	Lifetime limited warranty on frame 2-year limited warranty on electronics 18-months limited warranty on drive motors

Colors:



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

FEATURES

- 300 lbs. weight capacity
- Standard power elevating seat provides greater access to objects out of reach
- Powerful two-motor, patented mid-wheel drive design offers outstanding maneuverability
- Low center of gravity design delivers a very smooth and stable ride
- 12" knobby drive wheels and optimum weight distribution provide superior traction
- Speed range adjustable up to 4.5 mph
- Standard lap belt
- On-board battery charger

OPTIONS

- Contoured high-back seat
- Solid seat pan
- Solid tires
- Cup holder
- Heavy-duty dust cover
- Swing-away joystick
- Swing-away footrests
- Elevating leg rests
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Rear basket
- Angle-adjustable footplates
- Saddlebag
- Off-board charger

