



# Jazzy<sup>®</sup>

## SELECT

The easily maneuverable, Jazzy Select<sup>®</sup> is equipped with performance features like Active-Trac<sup>®</sup> Suspension and in-line motor technology for enhanced efficiency, torque and range for high performance. With features like a Contour high-back seat with a headrest, full length armrests and black, non-scuffing tires with silver rims, the Jazzy Select is a superb overall value.



**Pride<sup>®</sup>**  
Mobility Products Corp.  
— LIVE YOUR BEST<sup>®</sup> —

## SPECIFICATIONS

<b>Weight capacity</b>	300 lbs.
<b>Maximum speed</b>	Up to 4 mph
<b>Ground clearance</b>	1.875" at center frame
<b>Turning radius</b>	20"
<b>Overall length</b>	35.5" without front riggings
<b>Overall width</b>	23.75"
<b>Drive wheels</b>	10"
<b>Front casters</b>	5" solid
<b>Rear casters</b>	6" solid
<b>Suspension</b>	Active-Trac® Suspension
<b>Drivetrain</b>	Two-motor, in-line, Mid-Wheel 6®
<b>Braking system</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard electronics</b>	40 amp, PG GC 3 controller
<b>Specialty controls</b>	N/A
<b>Battery charger</b>	Off-board, 3 amp
<b>Per-charge range</b>	Up to 15 miles
<b>Battery requirements</b>	(2) 12 volt, deep cycle
<b>Battery size</b>	U-1
<b>Battery weight</b>	24.5 lbs. each
<b>Base weight</b>	99 lbs.
<b>Standard seat weight</b>	42.5 lbs. (Contour Seat)
<b>Max. seat dimensions</b>	20" x 20"
<b>Warranty</b>	5-year limited warranty on frame; 13-month limited warranty on electronics; 13-month limited warranty on drive motors

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

## MODELS

Jazzy Select 2-SS	(Code: K0822)
Jazzy Select 2-C	(Code: K0823)

## FEATURES

- New black, non-scuffing tires with silver rims
- Larger foot platform
- 40 amp, PG GC 3 controller
- Built with ease of service in mind
- Easy side access to batteries
- Toolless front panel access to electronics to minimize service time
- Active-Trac® Suspension
- Fold-flat, depth-adjustable high-back

## OPTIONS

- High-back seat with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

## RECOMMENDED LIFTS

**Silver Star®**  
**Lifts & Ramps**

- Backpacker® Plus
- Backpacker® MV
- Outlander
- Silver Boom 300



## COLOR CHOICES



Jazzy Red



Jazzy Blue

**PRIDE®**  
PERFORMANCE MOBILITY  
— LIVE YOUR BEST® —

(US) 800-800-8586 • (Canada) 888-570-1113 • [www.pridemobility.com](http://www.pridemobility.com)