

Jazzy[®]

SELECT 7

ULTRACOMPACT FRAME FOR
SMOOTH MANEUVERABILITY

Jazzy[®] Select[™] 7

The Jazzy[®] Select[™] 7 offers full-size power chair options and performance on an ultracompact frame for smooth maneuverability. A tight turning radius of only 20.5" and Active-Trac[®] Suspension give the Jazzy Select 7 strong outdoor performance without compromising superb indoor maneuverability.

#1 Physician Approved

JAZZY[®] SELECT[™] 7
shown in Jazzy Blue with high-
back seat, PG GC controller and
foot platform



Pride[®]
Mobility Products Corp.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

JZ SEL 7 SS - 08/07



The Jazzy® Select™ 7 series models include:

- Jazzy Select 7-2S-C
- Jazzy Select 7-2S-SS

Jazzy® Select™ 7

Weight capacity:	300 lbs.
Maximum speed:	Up to 4 mph
Ground clearance:	2.5"
Turning radius:	20.5"
Overall length:	36"
Overall width:	23"
Wheels:	Drive: 10"
Front anti-tips:	6" solid
Rear casters:	8" solid
Suspension:	Active-Trac®
Drive train:	Two-motor, mid-wheel drive
Braking system:	Intelligent braking (electronic, regenerative disc brakes)
Std. Electronics:	40 amp, PG GC controller
Specialty controls:	N/A
Battery charger:	On-board, 4 amp (standard)
Per-charge range:	Up to 15 miles
Battery requirements:	(2) 12 volt, deep cycle
Battery size:	U-1
Transport weight	
Battery weight:	24.5 lbs. each
Base weight:	77 lbs.
Seat weight:	37 lbs. (medium-back)
Max. Seat dimensions:	20" x 20"
Warranty:	5-year limited warranty on frame; 1-year limited warranty on electronics; 1-year limited warranty on drive motors

Colors:



FEATURES

- Standard Active-Trac® Suspension for varied terrain
- PG GC controller
- Powerful two-motor, mid-wheel drive design offers outstanding maneuverability
- Low center of gravity design delivers a very smooth and stable ride
- 10" knobby drive wheels and optimum weight distribution provide superior traction
- Speed range adjustable up to 4 mph
- On-board battery charger


OPTIONS

- Contoured high-back with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

Pride SilverStar Lifts & Ramps

- Backpacker® Plus
- Full Platform
- Power chair lift with straps
- Silver Boom 300



Backpacker® Plus



Full Platform

These lifts are compatible with the Jazzy® Select™ 7.