

Quantum® 1107

POWER CHAIRS

Quantum® 1107 shown with
Synergy® Manual Recline and
PG Remote Plus controller



Small frame. Big torque.

A highly stable and compact design gives the Quantum® 1107 outstanding rehab capability, while a powerful drivetrain delivers superb performance indoors and out. The Quantum 1107 is compatible with an impressive array of electronics and seating and positioning systems.



Specifications

DRIVE WHEELS	10"
REAR CASTERS:	6" solid
ANTI-TIPS:	5"
SUSPENSION	Advanced Active-Trac® with rear caster suspension
MAXIMUM SPEED ¹	Up to 5 mph
BRAKES	Intelligent Braking (electronic, regenerative disc brakes)
GROUND CLEARANCE	2.5"
TURNING RADIUS ³	21.875"
OVERALL SIZE ³	Length: 36.375" Width: 23.375"
SEATING	Synergy® Seat Adjustable Width: 10"-20" Adjustable Depth: 10"-20" TRU-Balance® 2 Power Positioning System ² Adjustable Width: 14"-20" Adjustable Depth: 14"-20" Specialty Seat
DRIVE TRAIN	2 motor, mid-wheel drive;
BATTERIES	Two 12 volt, deep cycle U-1
RANGE PER CHARGE ¹	Up to 15 miles
BATTERY CHARGER	5A, off-board
AVAILABLE ELECTRONICS	75A, Flight 70A, PG VSI 70A, Q-Logic
WEIGHT CAPACITY	300 lbs. maximum (user weight)
BASE WEIGHT	112 lbs.
BATTERY WEIGHT	24.5 lbs. (each)
WARRANTY	Lifetime limited warranty on frame 2-year warranty on electronics 18-month warranty on drive motor


Shroud Colors



Candy
Apple Red



Viper
Blue



Onyx
Black



Sunburst
Orange

(US) 866-800-2002
(Canada) 888-570-1113
www.quantumrehab.com



Features

- GLIDE motors for added power and torque
- Innovative battery boxes and front module eliminate battery and motor harnesses
- Advanced Active-Trac and rear caster suspension for enhanced performance on rough and uneven terrain, creating a smooth and stable ride
- Standard 5 amp, Lester Power Cube off-board charger

Options

- TRU-Balance® 2 Power Positioning System
- Q-Logic drive control system
- Contoured high-back seat
- Solid tires
- Cup holder
- Weather cover
- Swing-away joystick
- Elevating leg rests
- TRU-Comfort Seating
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Rear basket
- Angle-adjustable footplates
- 70° heavy-duty drop-in footrests

Code Coverage

Quantum 1107 2SP - SS | (Code: K0835)

FOOTNOTES:

- 1) Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) 200 lbs. weight capacity with TRU-Balance 2 Power Positioning System
- 3) Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %.