

QUANTUM®
FEEL THE POWER®
A Division of Pride Mobility Products® Corp.

Quantum® 1122

SERIES POWER CHAIRS

Quantum® 1122 shown with TRU-Balance® 2
Power Tilt & Recline and Remote Plus controller



The ultimate choice in mid-wheel drive.


The Quantum® 1122's unique design makes it a stable, stylish and high performing power chair. Its extensive adaptability and capability make The Quantum 1122 the ultimate choice in mid-wheel drive rehab power bases.



Specifications

DRIVE WHEELS	14"
REAR CASTERS:	8" solid
ANTI-TIPS:	6" solid, adjustable, front-mounted
SUSPENSION	Active-Trac® Suspension
MAXIMUM SPEED¹	Up to 6 mph (with high speed motors)
BRAKES	Intelligent Braking (electronic, regenerative disc brakes)
GROUND CLEARANCE	3.25"
TURNING RADIUS³	21"
OVERALL SIZE³	Length: 39.5" Width: 25.75"
SEATING	Synergy® Seat Adjustable Width: 10"-22" Adjustable Depth: 10"-20" TRU-Balance® 2 Power Positioning System Adjustable Width: 14"-22" Adjustable Depth: 14"-20" Contoured seating (max. size: 22" x 22") Manual recliner Solid seat pan option
DRIVE TRAIN	2 motor, mid-wheel drive; High-speed motor package (standard) High-torque / 4.5 mph (optional)
BATTERIES²	Two 12 volt, Group 24
RANGE PER CHARGE¹	Up to 20 miles
BATTERY CHARGER	8A, Off-board Lester Charger
AVAILABLE ELECTRONICS	70A, PG VSI 70A, PG Remote Plus
WEIGHT CAPACITY	300 lbs. maximum (user weight)
BASE WEIGHT	140 lbs.
BATTERY WEIGHT	53.5 lbs. (each)
WARRANTY	Lifetime limited warranty on frame 2-year warranty on electronics 18-month warranty on drive motor

Shroud Colors



Candy Apple Red Viper Blue Onyx Black Champagne Sunburst Orange

(US) 866-800-2002
(Canada) 888-570-1113
www.quantumrehab.com



Features

- 2-pole, patented two-motor, mid-wheel drive design provides superb power and acceleration
- Active-Trac® and rear caster suspension for outstanding outdoor performance
- Adjustable front anti-tips to meet your preferences and terrain challenges
- 8A, Lester off-board charger
- Easy front access to batteries

Options

- Solid tires
- Cup holder
- Heavy-duty dust cover
- Swing-away joystick
- Swing-away footrests
- Elevating leg rests
- Oxygen tank holder
- Cane holder
- Crutch holder
- Walker holder
- Rear basket
- Angle-adjustable footplates
- Saddlebag

Code Coverage

Q1122 3SP - SS	(Code: K0856)
Q1122 3SP - C	(Code: K0857)
Q1122 4MP - SS	(Code: K0884)

FOOTNOTES:

- 1) Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) NF-22 batteries are required with the power elevating seat option
- 3) Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %.