



# Quantum® 600 XL

SERIES POWER CHAIRS

Q600 XL shown with Synergy® Seating System and Flight controller



## Maximum stability.

The Quantum® 600 XL features our patent pending Mid-Wheel 6 design platform, an increased weight capacity of 400 lbs., and enhanced seating available up to 24" x 24". The 600 XL is the perfect heavy-duty choice for those who demand maximum stability and true all-around performance.



# Specifications

<b>DRIVE WHEELS</b>	14" solid
<b>OMNI-CASTERS</b>	<b>Front:</b> 5" solid <b>Rear:</b> 6" solid
<b>SUSPENSION</b>	ATX Suspension (Active-Trac <sup>®</sup> with extra stability)
<b>MAXIMUM SPEED <sup>1</sup></b>	Up to 4.5 mph
<b>BRAKES</b>	Intelligent Braking (electronic, regenerative disc brakes)
<b>GROUND CLEARANCE <sup>2</sup></b>	3"
<b>TURNING RADIUS <sup>2</sup></b>	20"
<b>OVERALL SIZE <sup>2</sup></b>	Length: 37.5" (without front riggings) Width: 23"
<b>SEATING</b>	Synergy <sup>®</sup> Seat Adjustable Width: 10"-24" Adjustable Depth: 10"-24" TRU-Balance <sup>®</sup> 2 Power Positioning System Adjustable Width: 14"-24" Adjustable Depth: 14"-24" High-back Seat (24" x 24" max. dim.)
<b>DRIVE TRAIN</b>	2-motor, GLIDE, Mid-Wheel 6 <sup>®</sup>
<b>BATTERIES</b>	Two 12 volt, NF-22
<b>RANGE PER CHARGE <sup>1</sup></b>	Up to 15 miles
<b>BATTERY CHARGER</b>	5A, off-board
<b>AVAILABLE ELECTRONICS</b>	75A Flight      70A PG Remote Plus 70A PG VSI
<b>WEIGHT CAPACITY</b>	400 lbs. maximum (user weight)
<b>BASE WEIGHT</b>	129 lbs.
<b>BATTERY WEIGHT</b>	38 lbs. (each)
<b>WARRANTY</b>	Lifetime limited warranty on frame 2-year warranty on electronics 18-month warranty on drive motor

# Shroud Colors

## METALLIC COLORS



## SOLID COLORS



# Features

- Mid-Wheel 6<sup>®</sup> allows six wheels on the ground for maximum stability
- OMNI-Casters (nylon, spherical-shaped casters) on front and rear prevent wheel hang-ups
- Side-mounted, easily accessible freewheel levers
- ATX Suspension (Active-Trac<sup>®</sup> with extra stability) incorporates front OMNI-Casters for enhanced performance over more varied terrain
- 5-amp, off-board charger
- 2-pole, GLIDE motors
- 400 lbs. weight capacity

# Options

- Swing-away joystick
- Multi-axis joystick mount
- Elevating leg rests
- Oxygen tank holder
- Neoprene leg rest sleeves
- Gel adductor buttons
- Angle-adjustable footplates
- Contoured seating
- High-mount foot platform
- Residual limb supports
- Contoured seating with solid seat pan

# Code Coverage

Q600 2HD - SS	(Code: K0824)
Q600 2HD - C	(Code: K0825)
Q600 2SPHD - SS	(Code: K0837)
Q600 2MPHD - SS	(Code: K0843)
Q600 3HD - SS	(Code: K0850)
Q600 3HD - C	(Code: K0851)
Q600 3SPHD - SS	(Code: K0858)
Q600 3SPHD - C	(Code: K0859)
Q600 3MPHD - SS	(Code: K0862)

## FOOTNOTES:

- 1) Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Due to manufacturing tolerances, this specification can be subject to variance of + or - 3%.



(US) 866-800-2002  
(Canada) 888-570-1113  
www.quantumrehab.com