

# Q610

## True All-Around Performance.

The compact Quantum® 610 is the ultimate choice for users who demand the superb stability of Mid-Wheel 6® Design and tight-quarter maneuverability. Highly adaptable, the Quantum 610 accepts a wide range of seating and electronics options while its powerful drivetrain delivers high-performance torque.



**Q610** shown with TRU-Balance® Power Positioning System and PG VR2 controller

## Specifications

<b>DRIVE WHEELS</b>	10" pneumatic or solid <sup>3</sup>
<b>OMNI-CASTERS</b>	<b>Front:</b> 5" solid <b>Rear:</b> 6" solid
<b>SUSPENSION</b>	ATX Suspension (Active-Trac® with extra stability)
<b>MAXIMUM SPEED<sup>1</sup></b>	Up to 4 mph
<b>BRAKES</b>	Intelligent Braking (electronic, regenerative disc brakes)
<b>GROUND CLEARANCE</b>	2.5"
<b>TURNING RADIUS<sup>2</sup></b>	21.25"
<b>OVERALL SIZE<sup>2</sup></b>	<b>Length:</b> 35.5" <b>Width:</b> 23.25"
<b>SEATING</b>	Synergy® Seat Adjustable Width: 10"-20" Adjustable Depth: 10"-20" TRU-Balance® Power Positioning System Adjustable Width: 14"-20" Adjustable Depth: 14"-20" TRU-Balance® Lift & Tilt Adjustable Width: 16"-20" Adjustable Depth: 16"-20" High-back Seat (20" x 20" max. dim.) Medium-back Seat (20" x 20" max. dim.)
<b>DRIVETRAIN</b>	2-motor, Mid-Wheel 6®
<b>BATTERIES</b>	Two 12 volt, deep cycle NF-22
<b>RANGE PER CHARGE<sup>1</sup></b>	Up to 20 miles
<b>BATTERY CHARGER</b>	5A, off-board
<b>AVAILABLE ELECTRONICS</b>	70A PG VR2 100A Q-Logic EX
<b>WEIGHT CAPACITY</b>	300 lbs. maximum (user weight) <sup>4</sup>
<b>BASE WEIGHT</b>	112 lbs.
<b>BATTERY WEIGHT</b>	24.5 lbs. (each)
<b>WARRANTY</b>	<b>K0822/K0823 Models:</b> Lifetime limited warranty on frame 13-month limited warranty on electronics 13-month limited warranty on drive motor  <b>K0835/K0836 Models:</b> Lifetime limited warranty on frame 2-year limited warranty on electronics 18-month limited warranty on drive motor

## Shroud Colors



## Features

- Mid-Wheel 6® allows six wheels on the ground for maximum stability
- Compatible with TRU-Balance® Power Positioning
- OMNI-Casters (nylon, spherical-shaped casters) on front and rear prevent wheel hang-ups
- Side-mounted, easily accessible freewheel levers
- ATX Suspension (Active-Trac® with extra stability) incorporates front OMNI-Casters for enhanced performance over more varied terrain
- Easy front access to batteries
- Powerful GLIDE motors for added power and torque

## Options

- TRU-Balance® Lift & Tilt
- TRU-Comfort Seating System
- Q-Logic EX Drive Controls
- LED full lighting package
- Power elevating seat
- Compact vent tray

## Code Coverage

Q610 2S - SS	(Code: K0822)
Q610 2S - C	(Code: K0823)
Q610 2SP - SS	(Code: K0835)
Q610 2SP - C	(Code: K0836)

### FOOTNOTES:

- 1) Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %.
- 3) Optional
- 4) 250 lbs. with TRU-Balance® Power Positioning Systems

**QUANTUM**  
**REHAB**®  
MOVING FORWARD

(US) 866-800-2002 ▪ (Canada) 888-570-1113 ▪ [www.quantumrehab.com](http://www.quantumrehab.com)