

VIBE®

REAR-WHEEL DRIVE POWER CHAIR

**NOW
WITH NEW
POWERFUL
GLIDE MOTORS**



ABOVE: Vibe shown with optional TRU-Balance® Cantilever seating system.
LEFT: Vibe shown with optional contoured high-back seat and swing-away footrests.

NEXT GENERATION REAR-WHEEL DRIVE THAT MANEUVERS LIKE A MID-WHEEL DRIVE...

800-800-8586 (US)
888-570-1113 (Canada)
www.pridemobility.com

Pride
Mobility Products Corp.

QUANTUM
REHAB
INNOVATIVE REHAB SOLUTIONS
A Division of Pride Mobility Products® Corp.

The Vibe® rear-wheel drive power chair goes far beyond conventional design to deliver unequalled versatility, striking style, and high-performance operation indoors and out. The first fully integrated mobility and positioning system, the Vibe's power base and seating seamlessly merge for optimal positioning and performance while its compact design allows it to maneuver like a mid-wheel drive power chair.

SPECIFICATIONS

Weight capacity:	300 lbs.
Maximum speed:	Up to 6 mph
Ground clearance:	3.75"
Turning radius:	23"
Overall length:	31.75"
Overall width:	25"
Wheels:	Drive: 14" pneumatic Casters: 9" pneumatic; 6", 8" and 9" solid Anti-tips: 4" solid
Suspension:	Sport-Trac
Drive train:	2-pole GLIDE, rear-wheel drive
Braking system:	Intelligent braking (electronic, regenerative disc brakes)
Std. Electronics:	50 amp, PG VSI controller
Specialty controls:	Yes
Battery charger:	Off-board, 8 amp (standard)
Per-charge range:	Up to 25 miles
Battery requirements:	(2) 12 volt deep cycle
Battery size:	Group 24 or Group 34
Battery weight:	Group 24: 53.5 lbs. Group 34: 43 lbs.
Base weight:	138 lbs.
Std. seat weight:	37 lbs. (Specialty Seat)
Max. seat dimensions:	20" x 20" (maximum rehab seating)
Warranty:	Lifetime limited warranty on frame 2-year limited warranty on electronics 18-months limited warranty on drive motors

COLOR CHOICES



Candy
Apple
Red



Viper
Blue



Champagne



Black

FEATURES

- High performance - high efficiency GLIDE motor package
- Easy operation results in a short learning curve for clients
- Narrow 25" width and compact 31.75" length
- Outstanding tight-space maneuverability
- Top speed up to 6 mph
- Sport-Trac Suspension
- 3-point seat mount for ease of adjustment and addition of optional components
- Power base and seating merged as a single unit for an optimal rehab package
- Easy rear access to batteries
- Group 34 low profile batteries with power positioning systems

OPTIONS

- Contoured medium/high-back seat
- TRU-Balance® Cantilever power elevating seat
- Synergy® Seat
- TRU-Balance® power positioning system
- Lap belt
- Flat-free tires
- Cup holder
- Heavy-duty dust cover
- 70° swing-away footrests
- Elevating leg rests
- Heel loops
- Oxygen tank holder
- Cane/crutch holder
- Rear basket
- Walker holder
- Swing-away joystick mount
- Saddlebag

(US) 800-800-8586
(Canada) 888-570-1113
www.pridemobility.com

Pride[®]
Mobility Products Corp.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.