



QUANTUM[®]
FEEL THE POWER[®]

A Division of Pride Mobility Products[®] Corp.

Quantum[®] Mid-Wheel Drive

REHAB POWER CHAIRS

Quantum[®] 1121



Exceptional Maneuverability

Quantum[®] Mid-Wheel Drive Power Chairs deliver exceptional maneuverability, high performance operation and superb rehab capability. Compatible with a wide range of rehab electronic and seating options including TRU-Balance[®] Power Positioning systems, Quantum Mid-Wheel Drive Power Chairs are designed to meet the diverse needs of the active rehab client.

Quantum[®]
1103 Ultra





1103 Ultra



1121

DRIVE WHEELS	10"	14"
REAR CASTERS:	8" solid	8" solid
ANTI-TIPS:	5" solid, adjustable, front-mounted	6" solid, adjustable, front-mounted
SUSPENSION	Active-Trac® Suspension	Active-Trac® Suspension
MAXIMUM SPEED ¹	Up to 4 mph	Up to 5 mph
BRAKES	Intelligent Braking (electronic, regenerative disc brakes)	Intelligent Braking (electronic, regenerative disc brakes)
GROUND CLEARANCE	2.5"	3.5"
TURNING RADIUS ²	20.75"	22.5"
OVERALL SIZE ²	Length: 35.5" Width: 23"	Length: 38.125" Width: 23.8"
SEATING	Synergy® Seat Adjustable Width: 10"-20" Adjustable Depth: 10"-20" TRU-Balance® Power Positioning System ³ Adjustable Width: 14"-20" Adjustable Depth: 14"-20" Contoured seating (max. size: 20" x 20") Manual recliner Solid seat pan option	Synergy® Seat Adjustable Width: 10"-20" Adjustable Depth: 10"-20" TRU-Balance® Power Positioning System Adjustable Width: 14"-20" Adjustable Depth: 14"-20" Contoured seating (max. size: 20" x 22") Manual recliner Solid seat pan option
DRIVE TRAIN	2 motor, mid-wheel drive	2 motor, mid-wheel drive
BATTERIES	Two 12 volt, U-1	Two 12 volt, NF-22
RANGE PER CHARGE ¹	Up to 15 miles	Up to 20 miles
STANDARD BATTERY CHARGER	5A, off-board	8A, off-board
AVAILABLE ELECTRONICS	50A, PG VSI controller 70A, PG Remote Plus controller	70A, PG VSI controller 70A, PG Remote Plus controller
WEIGHT CAPACITY	300 lbs. ³ maximum (user weight)	300 lbs. maximum (user weight)
BASE WEIGHT	92 lbs.	140 lbs.
BATTERY WEIGHT ⁴	24.5 lbs. (each)	38 lbs. (each)
WARRANTY	Lifetime limited warranty on frame 2-year warranty on electronics 18-month warranty on drive motor	Lifetime limited warranty on frame 2-year warranty on electronics 18-month warranty on drive motor
CODE COVERAGE	Q1103 Ultra 2SP - SS (Code: K0835) Q1103 Ultra 2SP - C (Code: K0836)	Q1121 3SP - SS (Code: K0856) Q1121 3SP - C (Code: K0857) Q1121 3MP - SS (Code: K0861)

Shroud Colors



Candy
Apple Red

Viper
Blue

Onyx
Black

Champagne

Sunburst
Orange



(US) 866-800-2002 • (Canada) 888-570-1113 • www.quantumrehab.com

- 1) Range and speed vary with user weight, terrain type, battery charge/condition and tire pressure.
- 2) Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %.
- 3) 200 lbs. with TRU-Balance power positioning systems.
- 4) Battery weight may vary based on manufacturer.