

JET™ 10

ULTRA



- Multiple footrest options
- Patented mid-wheel drive design provides outstanding maneuverability
- Speed range up to 4 mph
- On-board battery charger
- 10" pneumatic or flat-free drive wheels and optimum weight distribution provide superior traction
- Increased weight capacity of 300 lbs.

Jet™ 10 Ultra

Shown in Jet Blue with standard contoured medium-back gray vinyl seat with foot platform and non-programmable PG VSI controller.

The Jet™ 10 Ultra delivers a broad range of standard performance features like Active-Trac® Suspension and a powerful non-programmable 50-amp PG VSI controller on a very compact and maneuverable frame. Plus, the style is fresh, exciting and contemporary. These things combine with a truly economical price to make the Jet 10 Ultra an excellent overall value in personal mobility.



800-800-8586 (Exeter, PA) • 888-570-1113 (Canada) • www.pridemobility.com

FEATURES

- Compact size for superior handling in tight quarters
- Standard Active-Trac® Suspension
- 300 lbs. weight capacity
- Utilizes the trapeze bar seat mount for quick adjustability and easy interchange of a broad range of seats
- Convenient front battery access
- Curved rear caster beam provides excellent stability
- Maximum speed up to 4 mph
- Up to 25 miles of range per charge
- Standard non-programmable 50-amp PG VSI controller
- 18.5" turning radius

COLOR OPTIONS



Jet Red



Jet Blue

JET™ 10

ULTRA



Jet™ 10 Ultra

Shown in Jet Blue with standard contoured medium-back gray vinyl seat with foot platform and non-programmable PG VSI controller.

SPECIFICATIONS

Wheels:	10" solid or pneumatic
Rear Casters:	8" solid
Anti-tips:	5" solid, front-mounted
Suspension:	Active-Trac®
Maximum Speed:	Up to 4 mph
Brakes:	Intelligent braking (electronic, regenerative disc brakes)
Ground Clearance:	2.5"
Turning Radius:	18.5"
Overall Size:	Length: 35.5" Width: 23"
Seating Options:	Contoured medium-back, gray or beige vinyl (standard)
Drive Train:	Two-motor, mid-wheel drive
Battery:	Two 12 volt, U-I
Per-charge Range:	Up to 25 miles
Battery Charger:	4 amp, on-board (standard)
Electronics:	Non-programmable 50 amp, PG VSI
Weight Capacity:	300 lbs.
Weight:	Base: 97 lbs. Seat: 36.5 lbs.
Colors:	Jet Red and Jet Blue
Warranty:	5-year limited warranty on frame 1-year limited warranty on electronics 1-year limited warranty on drive motors

OPTIONS

- Contoured high-back, gray or beige vinyl seat
- Contoured seat with solid seat pan
- Lap belt
- Cup holder
- Heavy-duty dust cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Swing-away footrests
- Elevating leg rests
- Heel loops

800-800-8586 (Exeter, PA) • 888-570-1113 (Canada) • www.pridemobility.com

Pride
Mobility Products Corp.