

JET™ 2


Style, performance and value unite to make the Jet™ 2 one truly impressive power chair. Pride's patented JZMD® mid-wheel drive technology and sharp style bring the Jet 2 the performance and eye-catching look you demand. And the Jet 2's attractive price brings you the value you deserve.





	Jet™ 2	Jet™ 2 with HD package
Weight capacity:	300 lbs.	350 lbs.
Maximum speed:	Up to 4.5 mph	Up to 4.5 mph
Ground clearance:	3"	3"
Turning radius:	19.5"	19.5"
Overall length:	38"	38"
Overall width:	23.5"	24.5"
Wheels:	Drive: 14" pneumatic or solid	Drive: 14" pneumatic or solid
Front Anti-tips:	6" solid	6" solid
Rear Casters:	8" solid	8" solid
Suspension:	Limited	Limited
Drive train:	Two-motor, mid-wheel drive	Two-motor, mid-wheel drive
Braking system:	Intelligent braking (electronic, regenerative disc brakes)	Intelligent braking (electronic, regenerative disc brakes)
Std. Electronics:	Dynamic DL, 50 amp	PG VSI, 50 amp
Specialty controls:	N/A	N/A
Battery charger:	On-board, 4 amp	On-board, 4 amp
Per-charge range:	Up to 25 miles	Up to 25 miles
Battery requirements:	(2) 12 volt, deep cycle	(2) 12 volt, deep cycle
Battery size:	NF-22	NF-22
Transport weight		
Battery weight:	38 lbs. each	38 lbs. each
Base weight:	103 lbs.	104.5 lbs.
Std. seat weight:	37 lbs. (medium-back)	37 lbs. (medium-back)
Std. Seat dimensions:	18" x 18"	18" x 18"
Warranty:	5-year limited warranty on frame; 1-year limited warranty on electronics; 1-year limited warranty on drive motors	

Colors:



Jet Red



Jet Blue

FEATURES

- Available in Jet Red and Jet Blue
- Powerful two-motor, mid-wheel drive design
- Patented mid-wheel drive design offers outstanding maneuverability
- 14" knobby tires and optimum weight distribution provide superior traction
- Low center of gravity makes handling smooth and stable
- Speed range adjustable up to 4.5 mph
- Adjustable seat height and back angle
- On-board battery charger

OPTIONS

- Contoured high-back seat, gray or beige vinyl
- Lap belt
- Cup holder
- Heavy-duty dust cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Swing-away footrests
- Elevating leg rests
- Heel loops
- Height-adjustable armrests

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, weight, turning radius, and ground clearance specifications can be subject to variance of + or - 3 %.



- Full Platform
- Power chair lift with docking station
- Power chair lift with straps
- Backpacker



Backpacker



Full Platform

These lifts are compatible with the Jet™ 2 and Jet™ 2 HD.

Pride
Mobility Products Corp.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com