

# JET™ 2



## Jet™ 2

Shown in Jet Red with gray medium-back seat with foot platform and Dynamic controller.

Style, performance and value unite to make the Jet™ 2 one truly impressive power chair. Pride's patented JZMD® mid-wheel drive technology and sharp style bring the Jet 2 the performance and eye-catching look you demand. And the Jet 2's attractive price brings you the value you deserve.

**Pride**  
Mobility Products Corp.

800-800-8586 (Exeter, PA) • 888-570-1113 (Canada) • [www.pridemobility.com](http://www.pridemobility.com)

## FEATURES

# JET 2

- Available in Jet Red and Jet Blue
- Powerful two-motor, mid-wheel drive design
- Patented mid-wheel drive design offers outstanding maneuverability
- 14" knobby tires and optimum weight distribution provide superior traction
- Low center of gravity makes handling smooth and stable
- Speed range adjustable up to 4.5 mph
- Adjustable seat height and back angle
- On-board battery charger

## JET 2 WITH HD PACKAGE

- PG VSI, 50 amp controller (with standard 70 amp boost that engages automatically when needed)
- Available contoured seating up to 22" x 20"
- Up to 350 lbs. weight capacity
- 14" solid drive wheels (required)
- Base weight 104.5 lbs.
- Width 24.5"

### COLOR OPTIONS



## SPECIFICATIONS

<b>Drive Wheels:</b>	14" pneumatic or solid
<b>Rear Casters:</b>	8" solid
<b>Anti-tips:</b>	6" solid, front-mounted
<b>Suspension:</b>	Limited
<b>Maximum Speed:</b>	Up to 4.5 mph
<b>Brakes:</b>	Intelligent Braking (electronic, regenerative disc brakes)
<b>Ground Clearance:</b>	3"
<b>Turning Radius:</b>	19.5"
<b>Overall Size:</b>	Length: 38"    Width: 23.5"
<b>Seating Options:</b>	Contoured medium-back, gray or beige vinyl (standard)
<b>Drive Train:</b>	Two-motor, mid-wheel drive
<b>Battery:</b>	Two 12 volt, NF-22
<b>Per-charge Range:</b>	Up to 25 miles
<b>Battery Charger:</b>	4 amp, on-board
<b>Electronics:</b>	Dynamic DL, 50 amp
<b>Weight Capacity:</b>	300 lbs.
<b>Weight:</b>	Base: 103 lbs.    Seat: 37 lbs. (medium-back)
<b>Colors:</b>	Jet Red and Jet Blue
<b>Warranty:</b>	5-year limited warranty on frame 1-year limited warranty on electronics 1-year limited warranty on drive motors



## OPTIONS

- Contoured high-back seat, gray or beige vinyl
- Lap belt
- Cup holder
- Heavy-duty dust cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Swing-away footrests
- Elevating leg rests
- Heel loops
- Height-adjustable armrests

800-800-8586 (Exeter, PA) • 888-570-1113 (Canada) • [www.pridemobility.com](http://www.pridemobility.com)

**Pride**  
Mobility Products Corp.