

JET™ 7



Jet™ 7

Shown in Jet Red with contoured medium-back gray vinyl seat with foot platform and PG VSI controller.

The Jet™ 7 utilizes Pride's patented JZMD® mid-wheel drive design and a unique, easy-to-disassemble frame to bring you both outstanding indoor maneuverability and supreme portability. In addition to delivering a solid performance, the Jet 7 is priced to deliver something more...exceptional value.

Pride
Mobility Products Corp.

800-800-8586 (Exeter, PA) • 888-570-1113 (Canada) • www.pridemobility.com

FEATURES

- Frame disassembles for easy transport
- Available in Jet Red and Jet Blue
- Powerful two-motor, mid-wheel drive design
- Patented mid-wheel drive design provides outstanding maneuverability
- Low center of gravity creates a very smooth and stable ride
- 10" knobby drive wheels and optimum weight distribution provide superior traction
- Speed range adjustable up to 4 mph
- On-board battery charger

COLOR OPTIONS



JET™ 7



Jet™ 7

Shown in Jet Red with contoured medium-back gray vinyl seat with foot platform and PG VSI controller.

SPECIFICATIONS

Wheels:	10" pneumatic or solid
Rear Casters:	8" solid
Anti-tips:	6" solid; front-mounted
Suspension:	Limited
Maximum Speed:	Up to 4 mph
Brakes:	Intelligent Braking (electronic, regenerative disc brakes)
Ground Clearance:	2"
Turning Radius:	17"
Overall Size:	Length: 32.5" Width: 23"
Seating Options:	Contoured medium-back, gray or beige vinyl (standard)
Drive Train:	Two-motor, mid-wheel drive
Battery:	Two 12-volt, U-1
Per-charge Range:	Up to 25 miles
Battery Charger:	4 amp, on-board
Electronics:	50 amp, PG VSI
Weight Capacity:	300 lbs.
Weight:	Base: 89 lbs. Seat: 37 lbs. (medium-back)
Colors:	Jet Red and Jet Blue
Warranty:	5-year limited warranty on frame 1-year limited warranty on electronics 1-year limited warranty on drive motors

OPTIONS

- Lap belt
- Contoured high-back seat, gray or beige vinyl
- Cup holder
- Dust cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Height-adjustable armrests

800-800-8586 (Exeter, PA) • 888-570-1113 (Canada) • www.pridemobility.com

Pride
Mobility Products Corp.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.