



Pride® Lift Chairs
The Ultimate in Comfort & Style!

Pride®'s Classic Collection delivers superb value and excellent comfort with features like an innovative arm design and overstuffed fiber-filled backs. In addition, Classic Collection Lift Chairs feature patented, quiet and smooth operation and stylish designs to add to any home décor.

Classic COLLECTION

Formerly the
Casual Line

www.pridemobility.com

1-800-800-8586 (U.S.) 1-888-570-1113 (Canada)

Pride®
Mobility Products Corp.

CLASSIC COLLECTION MODELS



CL-10

2-Position
Medium to Large
Button-back Design
5'2" - 5'10"

CL-15

3-Position
Medium to Large
Button-back Design
5'2" - 5'10"

CL-20

2-Position
Medium to Large
Sewn Pillow-back Design
5'4" - 6'0"

CL-30

3-Position
Medium to Large
Chaise Lounger
Split-back Design
5'4" - 6'0"

SPECIFICATIONS

	CL-10	CL-15	CL-20	CL-30
Seat To Floor:	20"	20"	20"	20"
Seat Depth:	20"	20"	20"	20"
Seat Width:	20"	20"	20"	22"
Top of Back to Seat:	25"	27"	24.5"	26"
Distance From Wall:	12"	18"	11"	12"
Overall Chair Width:	33.5"	33.5"	32.5"	34"
Chair Weight	110 lbs.	119 lbs.	114 lbs.	118 lbs.

BENEFITS

- Patented, quiet & smooth lift system
- 375 lbs. weight capacity on CL-20 & CL-30; 325 lbs. weight capacity on CL-10 & CL-15
- Standard head and arm covers
- Two-year warranty
- Overstuffed blown fiber-filled back for comfort & support
- Integrated battery backup in case of power failure
- Heavy-duty lift actuator
- Fire retardant (California specifications; bulletin 117)
- Featuring Reflex™ environmentally friendly foam


FABRIC SELECTION




OPTIONAL FABRICS



PADDED-SUEDE



MICRO-SUEDE



LEVANTE VINYL



OPTIONS*

- *Head & lumbar pillows provide extra support & comfort*
- *Heat & massage for the ultimate relaxation experience (CL-20 & CL-30 only)*

* Additional charges apply for all options

Additional lead-time may apply for some options.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.