



Pride® Lift Chairs
The Ultimate in Comfort & Style!



Comfort and style complement one another perfectly in Pride's Heritage Collection Lift Chairs. With an array of plush fabrics to fit any décor and a superior design that provides exceptional strength, Heritage Collection chairs deliver style and durability that will last for years to come. Plus, Pride's Heritage Collection offers numerous standard features, superb comfort, and quiet and smooth performance, all at a great value.

HERITAGE COLLECTION

Formerly the
Gentle Line

www.pridemobility.com

1-800-800-8586 (U.S.) 1-888-570-1113 (Canada)

Pride®
Mobility Products Corp.

Pride® Lift Chairs

The Ultimate in Comfort & Style!

HERITAGE COLLECTION

HERITAGE COLLECTION MODELS



GL-58

2-Position
Medium to Large
5'2" - 5'10"

GL-358PW

3-Position
Petite Wide
Chaise Lounger
5'3" - Below

GL-358S

3-Position
Petite Chaise
Lounger
5'3" - Below

GL-358M

3-Position
Medium to Large
Chaise Lounger
5'4" - 6'0"

GL-358L

3-Position
Large to Extra-Large
Chaise Lounger
5'9" - 6'2"

GL-358XL

3-Position
Extra-Large
Chaise Lounger
5'6" - 6'1"

GL-358XXL

2-Position
Extra-Large
Chaise Lounger
5'6" - 6'1"

500 LBS. WEIGHT
CAPACITY

600 LBS. WEIGHT
CAPACITY

SPECIFICATIONS

	GL-58	GL-358PW	GL-358S	GL-358M	GL-358L	GL-358XL	GL-358XXL
Seat To Floor:	20"	19"	19"	20"	21"	20"	20"
Seat Depth:	20"	19"	18"	20"	22"	21"	22"
Seat Width:	20"	22"	18"	20"	22"	26"	30"
Top of Back to Seat:	24.5"	23.5"	23.5"	24.5"	25.5"	29"	29"
Distance From Wall:	11.5"	17"	17"	17"	17.5"	18"	18"
Overall Chair Width:	35"	39"	35"	35"	37"	39.5"	43"
Chair Weight	108 lbs.	118 lbs.	118 lbs.	117 lbs.	120 lbs.	159 lbs.	168 lbs.

FABRIC SELECTION



ALMOND



OATMEAL



MIDNIGHT
BLUE



EMERALD



MERLOT

OPTIONAL FABRICS



RASPBERRY



SAGE



MOCHA



GLACIER



ROSE TWIST

PADDED-SUEDE



BEACH



OLIVE



NEW KHAKI



NEW LATTE



BLUE



EXPRESSO



NAVY



KEYLIME

LEVANTE VINYL



BURGUNDY



ROYAL BLUE



MUSHROOM



NEW
CHESTNUT

BENEFITS

- Heavy-duty lift actuator with integrated battery backup
- 375 lbs. weight capacity (GL-358XL supports up to 500 lbs.; GL-358XXL supports up to 600 lbs.)
- Patented, quiet & smooth lift system
- Lifetime warranty on steel lift frame
- Seven-year limited warranty
- Head and arm covers standard
- Fire retardant (California specifications; bulletin 117)
- Single-switch hand control
- Standard dual pockets
- Featuring Reflex™ environmentally friendly foam

OPTIONS*

- *Massage & heat*
- *Head pillow*
- *Lumbar pillow*

Additional lead-time may apply for some options.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

www.pridemobility.com

1-800-800-8586 (U.S.) 1-888-570-1113 (Canada)

Pride®
Mobility Products Corp.