

# Stylus®

# LS

R E D E S I G N E D F O R 2 0 1 1

## Quality and Comfort... ...for Life.

The redesigned Stylus® LS meets the needs of active clients with a feature-rich and cost-effective package. All-new features include smooth drive tires and casters for increased ease of use, fluid-resistant, tension-adjustable seat upholstery and a lighter overall weight. The Stylus LS is a highly durable, highly adjustable manual wheelchair that delivers unmatched value.



**STYLUS® LS**  
shown with Synergy®  
Simplicity cushion

## SPECIFICATIONS

<b>PRODUCT WEIGHT</b>	32.5 lbs. without leg rests <sup>1</sup>
<b>WEIGHT CAPACITY</b>	250 lbs.
<b>FRAME MATERIAL</b>	Lightweight aluminum
<b>SEAT-TO-FLOOR HEIGHT<sup>2</sup></b>	17.5" with 24" wheels and 8" casters 16.5" with 22" wheels and 6" casters
<b>SEAT WIDTHS</b>	16", 18" and 20"
<b>SEAT DEPTHS</b>	16", 18" and 20" <sup>3</sup>
<b>BACK HEIGHTS</b>	17" - 19"
<b>REAR WHEELS</b>	22" or 24" composite mag wheels with low-profile, smooth polyurethane tires and aluminum handrims
<b>WHEEL LOCKS</b>	Push-to-lock
<b>CASTERS</b>	6" or 8" smooth polyurethane
<b>ANTI-TIP WHEELS</b>	Removable, adjustable
<b>WARRANTY</b>	Lifetime limited on side-frames and cross-frames; 1-year limited on materials

## RECOMMENDED LIFTS

- Caddy
- Silver Boom 300



**Silver Star**  
**Lifts & Ramps**

### FOOTNOTES:

- <sup>1</sup>) Standard 16" x 16" seat size and flip-up, desk-length arms
- <sup>2</sup>) Height varies from 16.5" - 19.5" based on wheel size
- <sup>3</sup>) Optional with extension kit

## FEATURES

- Fluid resistant, tension-adjustable seat upholstery
- Quick-release, removable rear wheels
- Low-profile, smooth 24" x 1" rear wheels
- Long handle wheel locks ensure comfort and security
- Industry-compatible leg rests with aluminum footplates for durability
- Quick-adjust footrests (no tools required)

## OPTIONS

- Synergy® Cushions
- Synergy® Shape Back
- Angle-adjustable back canes
- Tension-adjustable back upholstery
- Adjustable rear anti-tips
- Oxygen holder
- Elevating leg rests with calf pads
- Height-adjustable armrests, desk or full length
- Flip-up, desk length armrests
- Pull-to-lock wheel locks
- Wheel lock extensions
- Ultra-hemi package - 20" rear wheels and 5" casters with a 15" seat-to-floor height
- Lap belt
- 2" seat depth extension kit

**Pride**  
Mobility Products Corp.  
— LIVE YOUR BEST® —

(US) 800-800-8586 ▪ (Canada) 888-570-1113 ▪ [www.pridemobility.com](http://www.pridemobility.com)

The information contained herein is correct at the time of publication.  
We reserve the right to alter specifications without prior notice.