

	Revo® 3-wheel	Revo® 4-wheel
Model number:	SC63	SC64
Weight capacity:	300 lbs.	300 lbs.
Maximum speed <sup>1</sup> :	Up to 4.5 mph	Up to 4.5 mph
Ground clearance <sup>3</sup> :	2.25"	2.25"
Turning radius <sup>3</sup> :	38"	50"
Overall length <sup>3</sup> :	43.25"	44.75"
Overall width <sup>3</sup> :	21.5"	21.5"
Wheels		
Front:	2" x 8" solid	2" x 8" solid
Rear:	3" x 9" solid	3" x 9" solid
Anti-tips:	Rear	Rear
Suspension:	N/A	N/A
Transport weight		
Total w/o batteries <sup>3,6</sup> :	106.5 lbs.	125 lbs.
Heaviest piece when disassembled:	48 lbs. (rear section)	51.5 lbs. (front section)
Drivetrain:	Sealed transaxle, 24 volt DC motor	
Dual braking system:	Electronic, regenerative and electromechanical (electronic, regenerative disc brakes)	
Power elevating seat:	N/A	
Std. seating:	Foldable molded plastic (17"W x 17"D)	
Battery charger:	Off-board, 3 amp	
Per-charge range <sup>1,2</sup> :	Up to 20 miles	
Battery requirements:	(2) 12V deep cycle (AGM or gel cell type recommended)	
Battery & weight size <sup>3,4,5</sup> :	U-1 (24.5 lbs.each)	
Warranty:	3-year limited warranty	

Colors:



Red



Blue

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

<sup>2</sup> Tested in accordance with ASNI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

<sup>3</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

<sup>4</sup> AGM or Gel-Cell type required.

<sup>5</sup> Battery weight may vary based on battery manufacturer.

<sup>6</sup> Includes standard seat and weight.

## Benefits

- Exclusive black, non-scuffing tires
- Innovative interlocking drop-in battery packs for easy removal/installation
- Easily disassembles into five lightweight, easy-to-manage pieces
- 300 lbs. weight capacity
- State-of-the-art, microprocessor-based motor controller
- Seat offers height adjustment, 360° rotation and flip-up armrests for optimal user comfort
- Adjustable tiller for maximum comfort
- Headlight and taillight for added safety
- Easily-accessible charger port and circuit breaker

## Options

- Forearm crutch holder
- Lap belt
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Weather cover
- Saddlebag
- Quad cane holder
- Cane, clip attachment

Silver Star Lifts & Ramps 

- Backpacker® Plus
- Full Platform
- Silver Boom 300
- Backpacker® MV



Backpacker® Plus



Full Platform

These Silver Star® Lifts are compatible with the Revo®.

**Pride®**  
Mobility Products Corp.  
— LIVE YOUR BEST™ —

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

**PRIDE®**  
PERFORMANCE MOBILITY  
— LIVE YOUR BEST™ —

SLEEK STYLING  
AND SUPER PORTABILITY

## REVO®

The **Revo**® delivers full-sized scooter style, performance and features in a compact, lightweight design that easily disassembles for transport while remaining highly affordable.

