



Victory[®] XL

Outdoor Scooter



THE NEW VICTORY[®] XL Scooter delivers powerful outdoor performance and sleek styling, combined with truly innovative standard features — making it one of the finest scooter values available. The Victory XL offers a top speed up to 5.8 mph, full front and rear dynamic suspension, and comes standard with a convenient on-board battery charger as well as a full lighting package that includes a new, angle-adjustable lower headlight that is positioned for better pathway illumination.

Stylish
design
&
premium
performance
at an
excellent
value.



800-800-8586 (U.S.)
888-570-1113 (Canada)
www.pridemobility.com

Pride[®]
Mobility Products Corp.

Specifications

| | |
|------------------|---------------|
| Model number | SC2700 |
| Weight capacity | 400 lbs. |
| Maximum speed | Up to 5.8 mph |
| Ground clearance | 5" |
| Turning radius | 69" |
| Overall length | 55" |
| Overall width | 27" |

Tires

| | |
|--------------------------|---|
| Front | 4" x 13" pneumatic |
| Rear | 4" x 13" pneumatic |
| Anti-Tips | Rear |
| Standard seating | CRS deluxe reclining high-back with sliders & adjustable headrest |
| Power elevating seat | Option: Yes |
| Suspension | Front and rear |
| Drive train | Sealed transaxle, 24 volt DC motor |
| Dual braking system | Electronic, regenerative and electromechanical |
| Battery charger | On-board, 3 amp |
| Per-charge range | Up to 30 miles |
| Battery requirements | (2) 12V deep-cycle (AGM or gel cell) |
| Compatible battery sizes | NF-22, 40 AH, U-1 |

Weights:

| | |
|----------------------------------|---|
| Batteries | U-1: 24.5 lbs. each 40 AH: 32.5 lbs. each NF-22: 38 lbs. each |
| Total w/o batteries | 217 lbs. |
| Heaviest piece when disassembled | 87.5 lbs. (rear section) |

Colors



Candy
Apple Red



Viper Blue

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

Benefits

- Full directional lighting package, including an angle-adjustable lower headlight for pathway illumination
- High-back reclining seat with adjustable height and headrest for comfortable positioning
- Seat features sliders and swivels for additional comfort adjustments
- On-board battery charger—charge it wherever you go; up to a 30-mile range per battery charge
- Rugged four-wheel platform for stable outdoor performance
- Full dynamic suspension for a smooth ride, even on aggressive terrain
- Extensive range of options and accessories available
- Increased weight capacity of 400 lbs.
- 3-year limited warranty

Options

- Forearm crutch holder
- Lap belt
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Heavy-duty dust cover
- Saddlebag
- Quad cane holder
- Jumbo front basket
- Solid tires

