

VICTORY[®] 10

THE ULTIMATE 10" SCOOTER
WITH LUXURY FEATURES AT
AN ECONOMY PRICE.

VICTORY[®] 10

The sleek, sporty **Victory[®] 10** Scooter delivers high-performance operation, all-new features and feather-touch disassembly. With a full complement of advanced standard features like exclusive low-profile tires, backlit battery gauge and wraparound delta tiller, the Victory[®] 10 is a greater value than ever before.



Victory® 10	3-wheel	4-wheel
Model number:	SC610	SC710
Weight capacity:	400 lbs.	400 lbs.
Maximum speed ¹ :	Up to 5.25 mph	Up to 5.25 mph
Ground clearance ² :	3" at midframe	3" at midframe
Turning radius:	42"	54"
Overall length ² :	45.25"	47"
Overall width ² :	22.25"	22.25"
Range per charge ^{1,5} :	Up to 15.5 miles	Up to 15.5 miles
Tire size:	Front: 3.5x10 solid Rear: 3.5x10 solid	Front: 3.5x10 solid Rear: 3.5x10 solid
Total weight without batteries ^{1,2} :	131.5 lbs.	132.5 lbs.
Heaviest piece when disassembled ^{2,6} :	Front: 54.5 lbs.	Front: 61 lbs.
Standard seating:	Type: Foldable molded plastic Weight: 28 lbs. Material: Black Vinyl Dimensions: 18" width 17" depth	
Seat-to-deck range ² :	17.5" - 19"	
Drivetrain:	Rear-wheel drive, sealed transaxle, 24-volt DC motor	
Dual braking system:	Electronic, regenerative and electromechanical	
Battery requirements ³ :	Two 12-volt, deep cycle Size: U-1 Weight ⁴ : 24.5 lbs. each	
Battery charger:	3-amp, off-board	
Warranty:	3-year limited	

Colors:



Candy Apple Red



Viper Blue

¹Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

²Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³AGM or Gel-Cell type required.

⁴Battery weight may vary based on manufacturer.

⁵Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance test conducted at maximum weight capacity.

⁶Includes standard seat weight.

NOTE: All specifications subject to change without notice.

Features

- Exclusive one-hand feather-touch disassembly provides the ultimate quick transport system
- Auto-connecting front-to-rear harness
- Delta tiller with wraparound handles
- Pride®'s exclusive stylish, lightweight, non-scuffing, black, low-profile wheels
- Backlit battery gauge
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Lightweight seat with viscoelastic foam inserts for increased comfort and sliders
- Easy-grip tiller adjustment knob
- Easy access tiller-mounted charger port

Options

- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Weather cover
- Saddlebag

Silver Star Lifts & Ramps

- Backpacker® Plus
- Full Platform
- Backpacker® MV
- Silver Boom 300



Backpacker® Plus



Full Platform

These lifts are compatible with the Victory® 10.



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %.

