

VICTORY®
9 PS

FULLY LOADED WITH AN
IMPRESSIVE LIST OF FEATURES

VICTORY® 9 PS

The **Victory® 9 PS** delivers a long list of exciting features including a power elevating seat, low-profile tires, feather-touch disassembly, backlit battery gauge and wraparound delta tiller. With all of these luxury touches standard, along with high-performance operation and exhilarating style, the Victory 9 PS is a truly exceptional value.

*Standard
with a
power
elevating seat*




Pride®
Mobility Products Corp.
— LIVE YOUR BEST™ —

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

Victory® 9 PS 3-wheel

Model number:	SC609PS
Weight capacity:	300 lbs.
Maximum speed ¹ :	Up to 5.5 mph
Ground clearance ² :	2" at midframe
Overall length ² :	43.5"
Overall width ² :	23"
Turning radius:	41"
Range per charge ^{1,5} : <small>(Theoretical value based on test standards)</small>	13 miles
Tire size:	Front: 8" solid Rear: 9" solid
Total weight without batteries ^{1,2} :	131 lbs.
Heaviest piece when disassembled ² :	Front: 59.5 lbs.
Standard seating:	Type: Foldable molded plastic Weight: 28 lbs. Material: Black Vinyl Dimensions: 18" overall width 17" overall depth
Seat-to-deck range ² :	17.25" - 23.25"
Seat-to-ground range ² :	23.5" - 29.5"
Drivetrain:	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual braking system:	Regenerative, and electromechanical
Battery requirements ³ :	Two 12-volt, deep cycle Size: U-1 Weight: U-1:24.5 lbs. each ⁴
Battery charger:	3-amp, off-board

Colors:



¹Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

²Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³AGM or Gel-Cell type required.

⁴Battery weight may vary based on manufacturer.

⁵Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance test conducted at maximum weight capacity.

NOTE: All specifications subject to change without notice.

Features

- Standard power elevating seat
- Exclusive, one-hand feather-touch disassembly provides the ultimate quick transport system
- Pride's exclusive stylish, lightweight, low-profile, non-scuffing wheels
- Auto-connecting front-to-rear harness
- Delta tiller with wraparound handles
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Includes pinchless hinges and eight-position seat rotation
- 300 lbs. weight capacity

Options

- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Weather cover
- Saddlebag

Silver Star Lifts & Ramps

- Backpacker® Plus
- Full Platform
- Backpacker® MV

