



LASER[®]

3 and 4-wheel Sport Scooters[®]



In addition to Turbo Mode and patent-pending Speed Control Technology (SCT), the Laser[®] is fully loaded with innovative standard features like an infinitely adjustable tiller, dynamic front and rear suspension, full lighting package, high-density front and rear bumpers, and our sleek, contoured Sport Scooter[®] high-back seat. All of these combine with beautiful styling to offer you a Sport Scooter that keeps up with your active life and complements your personal style.

Featuring
Turbo Mode
for reaching
speeds of
8 mph* in
straightaways.



Laser 4-wheel Sport Scooter shown in Viper Blue

* Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

800-800-8586 (U.S.)
888-570-1113 (Canada)
www.pridemobility.com

sport scooter *n.* (sport - skoot r): 1. A scooter that can go 8 mph in straightaways and due to Speed Control Technology, its speed is tightly controlled when turning, going down declines and laterally banking. 2. See Pride Mobility Products Corp.

Pride[®]
Mobility Products Corp.

Specifications

Model number(s)	3-wheel: SPSC300 4-wheel: SPSC340
Weight capacity	300 lbs.
Maximum speed	Up to 8 mph (in Turbo Mode)
Ground clearance	3.5"
Turning radius	3-wheel: 41" 4-wheel: 58"
Overall length	3-wheel: 48.5" 4-wheel: 50.25"
Overall width	26"

Tires

Front	3" x 10" pneumatic
Rear	4" x 10" pneumatic
Anti-tips	Rear
Standard seating	Sport Scooter high-back with headrest and sliders
Suspension	Front and rear
Drive train	Sealed transaxle, 24 volt DC motor
Dual braking system	Electronic, regenerative and electromechanical
Battery charger	On-board, 3 amp
Per-charge range	Up to 30 miles
Battery requirements	(2) 12V deep-cycle
Compatible battery sizes	32AH (U-1); 40AH

Weight

Battery	32AH: 24.5 lbs. 40AH: 32.5 lbs.
Total w/o batteries	3-wheel: 152.5 lbs. 4-wheel: 169.5 lbs.
Heaviest piece	63 lbs. (rear section)

Colors



Viper
Blue



Candy
Apple
Red



Black
Onyx

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

Benefits

- Turbo speed mode for up to 8 mph in straightaways
- Speed Control Technology (SCT) automatically reduces the speed during cornering, steep descents and lateral banking
- Dynamic front and rear suspension
- Full lighting package, including directional signals
- Easily accessible circuit breaker/charger port and freewheel lever
- Tiller infinitely adjusts to your most comfortable driving position
- Contoured, supportive high-back sport seat with headrest
- User-adjustable sliding seat

Options

- Power elevating seat
- Forearm crutch holder
- Lap belt
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Heavy-duty dust cover
- Saddle bag
- External battery charger
- Quad cane holder