

# Jazzy<sup>®</sup>

1107

# Jazzy<sup>®</sup>

1107



Jazzy<sup>®</sup> 1107  
Shown in Blue with  
medium-back seat,  
footboard and VSI controller



More Advanced Features, Excellent Performance & Super Portability. The Jazzy<sup>®</sup> 1107 brings the most advanced achievements of our engineering efforts together in one compact yet superbly balanced package that conveniently disassembles for transport. Direct connect components eliminate wiring harnesses and provide clean style, unquestionable convenience, and quick and easy disassembly.

Jazzy<sup>®</sup> 1107  
Shown in Red with  
medium-back seat,  
footboard and VSI controller



# Specifications:

<b>WEIGHT CAPACITY:</b>	136kg (21 STONE)
<b>MAXIMUM SPEED:</b>	Up to 6 km/h (4 mph)
<b>GROUND CLEARANCE:</b>	6cm (2.5")
<b>TURNING RADIUS:</b>	52cm (20.5") <i>Without Foot Riggings</i>
<b>LENGTH:</b>	93cm (36.5") <i>Without Foot Riggings</i>
<b>WIDTH:</b>	58cm (23")
<b>DRIVE WHEELS:</b>	25.5cm (10") pneumatic
<b>CASTORS:</b>	15cm (6") solid
<b>ANTI-TIPS:</b>	13cm (5") solid
<b>SUSPENSION:</b>	Advanced Active-Trac* with independent rear castors
<b>DRIVE TRAIN:</b>	2-motor, mid-wheel drive
<b>BRAKING SYSTEM:</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>ELECTRONICS:</b>	VSI controller
<b>STANDARD SEATING:</b>	Medium-back grey vinyl
<b>STANDARD FOOTREST:</b>	Footboard
<b>BATTERY CHARGER</b>	Off-board, 4 amp (standard)
<b>PER-CHARGE RANGE:**</b>	Up to 19km (12 miles) (max.)
<b>BATTERIES:</b>	Two 12 volt 30 amp
<b>BATTERY WEIGHT:</b>	11kg (24.5 lbs.) each
<b>BASE WEIGHT:</b>	51kg (112 lbs.)
<b>SEAT WEIGHT:</b>	27kg (41 lbs.) (medium-back)
<b>MAX. SEATING</b>	See prescription form
<b>DIMENSIONS:</b>	for seat sizing
<b>WARRANTY:</b>	Lifetime warranty on frame 1-year limited warranty on electronics 2-year limited warranty on drive train 18-month PG Drive control systems



- Locking integration of four components provides rigid frame stability



- Direct connectors eliminate wiring harnesses



- Quickly and easily disassembles for travel

# Features:

# Jazzy


1107

- Advanced Active-Trac® and rear castor suspension for a smoother ride over challenging terrain
- Quickly and easily disassembles for transportation
- Aluminum frame rails and centre locking mechanism provide a very rigid platform for ultimate stability
- Innovative battery boxes and front module eliminate battery and motor harnesses
- VSI controller
- Flip-up armrests

# Options:

- Solid tyres
- Cup holder
- Cane/crutch holder
- Rear basket
- Swing-away joystick
- Swing-away footrests
- Elevating leg rests
- Oxygen tank holder
- 80° heavy-duty, drop-in footrests

# Colours:



\*The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice.\*

\*\*Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.

\*\*\*Exact product configurations shown on this literature may not be available in all markets

- Active-Trac® & rear castor suspension



- Cup holder



- Cane/crutch holder

