

Jazzy[®]

SELECT 6

The Jazzy Select[®] 6 provides maximum front stability, combined with the advantages of in-line motor technology. Active-Trac[®] Suspension and 6" front and rear castors offer superior performance indoors and out. An impressive selection of standard convenience features make it the easiest power chair to use.



Pride[®]
Mobility Products Ltd.
— LIVE YOUR BEST —

SPECIFICATIONS

Weight capacity	21 Stone 6 lbs. (136 kg)
Maximum speed¹	Up to 4 mph (6.44 km/h)
Ground clearance	5.08 cm (2.0")
Turning radius²	55.56 cm (21.8750")
Overall length²	89.85 cm (35.3750") ⁴
Overall width²	59.69 cm (23.5")
Wheels	Drive: 25.4 cm (10") solid
Anti-tips	N/A
Castors	15.24 cm (6") solid (front and rear)
Suspension	Active-Trac [®]
Drivetrain	Two-motor, in-line, mid-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	40 amp, PG GC/ 60 amp VR2 (with power elevating seat)
Specialty controls	No
Battery charger	Off-board, 5 amp (standard)
Battery requirements	(2) 12 volt deep cycle
Battery size	(U-1) 30 AH
Battery weight³	11.11 kg (24.5 lbs.) each
Base weight	45.36 kg (100 lbs.)
Standard seat weight	16.78 kg (37 lbs.) high back
Max. seat dimensions	46 cm x 46 cm (18" x 18")
Warranty	Lifetime warranty on frame; 1-year limited warranty on electronics; 18-month warranty on PG Drive control system; 2-year limited warranty on drivetrain

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure.
- 2) Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.
- 3) Battery weight may vary based on manufacturer.
- 4) Without foot riggings.

FEATURES

- 15.24 cm (6") castors on front and rear for maximum stability
- In-line motor technology for increased efficiency, torque and range
- 40 amp, PG GC controller/ 60 amp VR2 (some models)
- Articulating rear castor beam
- Battery removal through side panel
- Easy access charger port and circuit breaker
- Service friendly access to controller
- Height and depth-adjustable foot platform

OPTIONS

- Power elevating seat
- Cup holder
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

COLOURS



Red



Blue



PRIDE[®]
PERFORMANCE MOBILITY
— LIVE YOUR BEST[®] —

32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL
(Tel) 01869 324 600 • (Fax) 01869 323 070 • www.pridemobility.co.uk