

## Pride® Lift Chairs

*Experience the Luxury*

# LC-101 Lift Chair



### A 3-Position Lift Chair that Delivers Superb Value and Excellent Comfort

The LC-101 from Pride® delivers superb value, excellent comfort and exceptional quality. The LC-101 features quiet and smooth operation and stylish designs to add to any home décor.

#### FABRIC SELECTION



Cocoa



Sandal

Note: Measurements may vary +/- 2.5 cm (1").

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

#### SPECIFICATIONS

##### LC-101

<b>Back Width:</b>	72 cm (28.5")
<b>Seat-to-Floor:</b>	47 cm (18.5")
<b>Seat Depth:</b>	51 cm (20")
<b>Seat Width:</b>	51 cm (20")
<b>Top of Back to Seat:</b>	70 cm (27.5")
<b>Weight Capacity:</b>	147 kg (325 lbs.)

*Experience the Luxury*



**Sawa Medical Supplies, s.a.r.l.**

Tel: 01/804.881 03/639.650

Raouche, Beirut Lebanon

Website: [www.sawamed.com](http://www.sawamed.com)

EEmail: [info@sawamed.com](mailto:info@sawamed.com)

**sawa**  **سَوَا** شركة  
medical supplies s.a.r.l. للمعدات الطبية ش.م.ل.