


Wrangler® 4-wheel

Model number:	PMV600
Weight capacity:	400 lbs.
Maximum speed ¹ :	Up to 10 mph
Ground clearance ² :	7"
Turning radius:	77"
Overall length ² :	63.5"
Overall width ² :	31"
Wheels Front: Rear: Anti-tips:	6" x 13" pneumatic 6" x 13" pneumatic N/A
Suspension:	Rear
Total weight w/o batteries ^{2,5} :	216 lbs.
Heaviest piece when disassembled ² :	No frame disassembly
Drive train:	Two in-line gear motors
Dual braking system:	Electronic, regenerative and electromechanical (electronic, regenerative disc brakes)
Power elevating seat:	Yes (option)
Std. Seating:	Deluxe contoured high-back with headrest
Battery charger:	On-board, 5 amp
Per-charge range ¹ :	Up to 20 miles
Battery requirements ³ :	(2) 12V deep cycle (AGM or gel cell type recommended)
Battery size & weight ⁴ :	Group 27; 63.5 lbs. each
Warranty:	3-year limited warranty

Colors:



Candy Apple Red



Viper Blue

Wrangler Benefits

- 200-amp, Curtis microprocessor-based PMV controller system provides optimal levels of torque and power
- Two rugged, independent in-line gear motors provide positive traction
- Variable speed up to 10 mph

Options

- Rear basket
- Safety flag
- Heavy-duty dust cover
- Weather breaker canopy
- Cup holder
- Rear basket

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

² Tested in accordance with ASNI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

⁴ AGM or Gel-Cell type required. Battery weight may vary due to manufacturer.

⁵ Includes standard seat & weight.

Wrangler®

STYLISH DESIGN &
PREMIUM PERFORMANCE

WRANGLER®

The **Wrangler®** is a powerful blend of rugged, heavy-duty performance with refined luxurious features. Whether you find yourself limited by a mobility impairment that threatens to interfere with your outdoor activities or you simply want a stylish, practical and durable mobility vehicle that is easy to use, quiet and comfortable, the Wrangler is the perfect choice.



Pride®
Mobility Products Corp.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com