

Backpacker[®]

MV



Now you can take your mobility on the road easily in your Sport Utility Vehicle with the Backpacker MV from Pride Silver Star[®] Lifts and Ramps. It's the ultimate solution to lifting your powered mobility device into your vehicle, engineered specifically to be a perfect fit for a wide variety of popular vehicles including SUVs. Like all Silver Star Lifts from Pride, the Backpacker MV features all-steel construction and a protective finish, along with an industrial rated 12 volt motor and battery pack, and easy-to-use controls to combine rugged durability with safe, simple and reliable operation.

Other features include:

- One piece, fully assembled drive-on/drive-off platform
- Removable, front drive panel for easier serviceability
- New, adjustable wheel chocks on the platform surface
- Single piece horizontal drive rail: Minimum welding points for noise reduction, increased rigidity, and more fluid movement
- New, angled, low-profile platform
- Easy-to-use 2 button hand control
- Can lift a wide range of Pride[®] Scooters and Power Chairs
- 159 kg (350 lbs.) lifting capacity
- 3-year transferable warranty
- Fits into SUVs with a 79 cm (31") opening height
- Manual override
- Improved micro switch protection
- False floor available for specified vehicles through special order to preserve vehicle aesthetics

Pride[™] Lifts & Ramps

20-24 Apollo Drive, Hallam, Victoria 3803

ACN 088 609 661 ABN 66 088 609 6611

(Tel) 03 8770 9600 • (Fax) 03 9703 2720 • www.pridemobility.com.au

Backpacker® MV Interior Lift



Backpacker® MV in lowest position.



Backpacker® MV fully loaded.

Total lift system weight	87 kg (191 lbs.) ¹
Scooter / power chair overall weight	159 kg (350 lbs.)
Overall platform length	112 cm (44") (104 cm (41") usable surface)
Maximum product width	68.5 cm (27")
Maximum product wheelbase	104 cm (41")
Maximum travel distance from the ground	94 cm (37")
Motor	Sealed 12 volt DC
Warranty	3-year transferable
Optional battery pack /integrated charging system	12 volt, on-board
Optional charger	110 volt

Can easily lift the following current Pride® Power Chairs²

- Jazzy® 600
- Jazzy® 600 XL
- Jazzy® 1113 ATS
- Jazzy® Select™
- Jazzy® Select™ 6
- Jazzy® Select™ 6 Ultra
- Jazzy® Select™ 7
- Jazzy® Select™ 14
- Jazzy® Select™ 14XL
- Jazzy® Select™ Traveller
- Jet™ 3 Ultra

Can easily lift the following current Quantum® Power Chairs²

- Quantum® 600
- Quantum® 600 XL
- Quantum® 610
- Quantum® 1103 Ultra
- Quantum® 1121
- Quantum® 6000
- Quantum® 6000 XL
- Quantum® 6000Z
- Quantum® R4000

Can easily lift the following current Pride® Scooters²

- Celebrity® X 3 & 4
- Go-Go® Ultra X 3 & 4
- Go-Go Elite Traveller® 3 & 4
- Go-Go Elite Traveller® Plus 3 & 4
- Go-Go Elite Traveller® Plus HD 3 & 4
- Legend® 3 & 4
- Maxima 3 & 4
- Rally® 3 & 4
- Revo® 3 & 4
- Sundancer®
- Victory® 10 3 & 4

- | |
|---|
| <p>¹ Weight does not include any accessories</p> <p>² Vehicle dependent</p> <p>³ False floor kits available for specified vehicles through special order</p> |
|---|

Note: The **Backpacker® MV** can lift all of the listed power chairs with contoured, Synergy® or Specialty seating. Please call Customer Service when requesting information about lifting with a power positioning system. Can fit into most SUV's and mini-vans with little or no vehicle modification, full-size vans require modification.

Pride™ Lifts & Ramps

20-24 Apollo Drive, Hallam, Victoria 3803

ACN 088 609 661 ABN 66 088 609 6611

(Tel) 03 8770 9600 • (Fax) 03 9703 2720 • www.pridemobility.com.au



Quality System
Quality Endorsed Company
ISO 9001 OEC20029
SAI Global