

Jazzy[®]

SELECT SERIES

Jazzy Select[®] Power Chairs offer the advantages of inline motor technology for enhanced efficiency, torque and range. The highly manoeuvrable Jazzy Select is loaded with performance and convenience features, while the Jazzy Select 6 provides maximum front stability for superior performance indoors and out. The Jazzy Select 14 XL completes the package, with a powerful 60 amp PG VR2 controller.

Jazzy Select[®] 6
shown in Jazzy Red



Jazzy Select[®]
shown in Jazzy Red



Jazzy Select[®] 14 XL
shown in Jazzy Red



21 Healey Rd., Dandenong, Victoria 3175

ACN 088 609 661 ABN 66 088 609 661

(Tel) 03 9706 4611 • (Fax) 03 9706 4622 • www.pridemobility.com

Pride[®]
Mobility Products
Australia Pty. Ltd.

SPECIFICATIONS	Jazzy Select®	Jazzy Select® 14XL	Jazzy Select® 6
Weight capacity	136 kg (300 lbs.)	181 kg (400 lbs.)	136 kg (300 lbs.)
Maximum speed	Up to 6 km/h (4 mph)	Up to 6.5 km/h (4mph)	Up to 6 km/h (4 mph)
Ground clearance	4.76 cm (1.875")	7.6 cm (3")	5 cm (2")
Turning radius	51 cm (20")	53 cm (21")	52 cm (20.5")
Overall length¹	90 cm (35.5")	96.5 cm (38")	93 cm (36.625")
Overall width	60 cm (23.75")	63 cm (24.75")	60 cm (23.5")
Wheels	Drive: 25 cm (10")	Drive: 35.5 cm (14")	Drive: 25 cm (10") solid
Front castors	13 cm (5") solid anti-tips	15 cm (6") solid anti-tips	15 cm (6") solid
Rear castors	15 cm (6") solid	20 cm (8") solid	15 cm (6") solid
Suspension	Active-Trac®	Limited	Active-Trac®
Drive train	Two-motor, in-line, mid-wheel drive	Two-motor, mid-wheel drive	Two-motor, in-line, mid-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)		
Standard electronics	40 amp, PG GC 2	60 amp, PG VR2	40 amp, PG GC 2
Specialty controls	N/A		
Battery charger	Off-board, 3 amp	Off-board, 4-amp	Off-board, 3 amp
Per-charge range	Up to 24 km (15 miles)		
Battery requirements	(2) 12 volt, deep cycle		
Battery size	35 AH	35 AH	35 AH
Battery weight	11 kg (24.5 lbs.) each	11 kg (24.5 lbs.) each	11 kg (24.5 lbs.) each
Base weight	45 kg (99 lbs.)	48 kg (106.5 lbs.)	45 kg (100 lbs.)
Standard seat weight	19 kg (42.5 lbs.) (high-back)	19 kg (42.5 lbs.) (high-back)	19 kg (42.5 lbs.) (high-back)
Max. seat dimensions	46 cm x 46 cm (18" x 18")	46 cm x 51 cm (18" x 20")	46 cm x 51 cm (18" x 20")
Warranty	2-year limited		

COLOUR CHOICE



Jazzy Red

1) Without front riggings



FEATURES (Jazzy Select®)

- Dual in-line motors for enhanced efficiency, torque, range and performance
- New 40 amp, PG GC 2 controller
- Build with ease of service in mind
- Easy side access to batteries
- Tool-less front panel access to electronics to minimize service time

FEATURES (Jazzy Select® 14 XL)

- Patented mid-wheel drive design for tight space manoeuvrability
- 35.5 cm (14") knobby tyres and optimum weight distribution provide superior traction
- Upgradeable to 55 AH batteries for maximum range per charge
- New 60 amp, PG VR2 controller
- Off-board, 4 amp charger

FEATURES (Jazzy Select® 6)

- 15 cm (6") casters on front and rear for maximum stability
- In-line motor technology for increased efficiency, torque and range
- Articulating rear caster beam
- Easy access charger port and circuit breaker
- Height and depth-adjustable foot platform
- Built with ease of service in mind

OPTIONS

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder

21 Healey Rd., Dandenong, Victoria 3175

ACN 088 609 661 ABN 66 088 609 661

(Tel) 03 9706 4611 • (Fax) 03 9706 4622 • www.pridemobility.com

The information contained herein is correct at the time of publication; due to the Pride's commitment to constant improvement and development, We reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

PRIDE[®]
PERFORMANCE MOBILITY
— LIVE YOUR BEST™ —