

Pride® Lift Chairs
The Ultimate in Comfort & Style



Shown in Sandal



Shown in Fern

LC-101

Lift Chair


A 3-Position Lift Chair with classic style and exceptional quality and performance

Pride's LC-101 delivers superb value, excellent quality and luxurious comfort. The LC-101 features quiet and smooth operation and stylish design to add to any home décor.

SPECIFICATIONS

	LC-101
Seat to Floor:	47 cm (18.5")
Seat Depth:	51 cm (20")
Seat Width:	51 cm (20")
Top of Back to Seat:	70 cm (27.5")
Distance From Wall:	45.75 cm (18")
Overall Chair Width:	83 cm (32.75")
Weight Capacity:	147 kg (325 lbs.)
Available Back Styles:	Button-back

FABRIC SELECTION (MICRO-SUEDE)



Cocoa



Fern



Sandal

Note: Measurements may vary +/- 2.5 cm (1").

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

Pride®
Mobility Products
Australia Pty. Ltd.

21 Healey Rd.
Dandenong, Victoria 3175
www.pridemobility.com
ACN 088 609 661 ABN 66 088 609 661

Tel: 03 9706 4611
Fax: 03 9706 4622

Experience the Luxury

