

QUANTUM®
FEEL THE POWER®

Quantum® 1122

SERIES POWER CHAIRS

Quantum® 1122 shown with TRU-Balance® 2
Power Tilt & Recline and Remote Plus controller



The ultimate choice in mid-wheel drive.

The Quantum® 1122's unique design makes it a stable, stylish and high performing power chair. Its extensive adaptability and capability make The Quantum 1122 the ultimate choice in mid-wheel drive rehab power bases.



Pride®
Mobility Products
Australia Pty. Ltd.

Specifications

DRIVE WHEELS	35.5 cm (14")
REAR CASTORS:	20 cm (8") solid
ANTI-TIPS:	15 cm (6") solid, adjustable, front-mounted
SUSPENSION	Active-Trac® Suspension
MAXIMUM SPEED ¹	Up to 9.5 km/h (6 mph) (with high speed motors)
BRAKES	Intelligent Braking (electronic, regenerative disc brakes)
GROUND CLEARANCE	8.25 cm (3.25")
TURNING RADIUS ³	53 cm (21")
OVERALL SIZE ³	Length: 100 cm (39.5") Width: 65 cm (25.75")
SEATING	Synergy® Seat Adjustable Width: 25 cm - 61 cm (10"-24") Adjustable Depth: 25 cm - 61 cm (10"-24") TRU-Balance® 2 Power Positioning System Adjustable Width: 35.5 cm - 61 cm (14"-24") Adjustable Depth: 35.5 cm - 61 cm (14"-24") Contoured seating (max. size: 61 cm x 61 cm (24" x 24")) Manual recliner Solid seat pan option
DRIVE TRAIN	2 motor, mid-wheel drive; High-speed motor package High-torque / 7.25 km/h (4.5 mph) (optional)
BATTERIES ²	Two 12 volt, Group 24
RANGE PER CHARGE ¹	Up to 32 km (20 miles)
BATTERY CHARGER	8A, Off-board charger
AVAILABLE ELECTRONICS	70A, PG VSI 70A, PG Remote Plus
WEIGHT CAPACITY	136 kg (300 lbs.)/ XL Package: 181 kg (400 lbs.) maximum (user weight)
BASE WEIGHT	63.5 kg (140 lbs.)
BATTERY WEIGHT	24 kg (53.5 lbs.) (each)
WARRANTY	2-year limited

Features

- 2-pole, patented two-motor, mid-wheel drive design provides superb power and acceleration
- Active-Trac® and rear castor suspension for outstanding outdoor performance
- Adjustable front anti-tips to meet your preferences and terrain challenges
- 8A, off-board charger
- Easy front access to batteries

Options

- Solid tyres
- Cup holder
- Heavy-duty dust cover
- Swing-away joystick
- Swing-away footrests
- Elevating leg rests
- Oxygen tank holder
- Cane holder
- Crutch holder
- Walker holder
- Angle-adjustable footplates

Shroud Colours



Candy
Apple Red

Viper
Blue

Onyx
Black

Champagne

Sunburst
Orange

FOOTNOTES:

- 1) Range and speed vary with user weight, terrain type, battery charge, battery condition and tyre pressure.
- 2) NF-22 batteries are required with the power elevating seat option
- 3) Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %.

"The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice".

"Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement".