

QUANTUM

for kids



The Q6 Power Chair Series and the Dynamo® 1107 deliver grown up power and performance that's up for any challenges a kid can put it through. Limitless rehab capability, powerful motors and the unmatched manoeuvrability of our patent pending Mid-Wheel 6® design make Q6 Series power chairs an excellent choice for the classroom, playground and beyond. All Q6 Series power chairs and the Dynamo 1107 offer exceptional versatility and adjustability, allowing them to grow with kids to meet their changing needs.

Q610



Q600



Q6000Z



Dynamo®
1107



THE POWER YOU TRUST

See the Q6 Series in action at:
www.quantumrehab.com




	Q610	Q600	Q6000Z	Dynamo® 1107
Drive wheels	25.4 cm (10")	35.5 cm (14")	35.5 cm (14")	25.4 cm (10")
Front wheels	13 cm (5") solid castors	13 cm (5") solid castors	15 cm (6") solid castors	13 cm (5") solid anti-tips
Rear wheels	15 cm (6") solid castors	15 cm (6") solid castors	15 cm (6") solid castors	15 cm (6") solid castors
Suspension	Active-Trac® ATX Suspension	Active-Trac® ATX Suspension	Active-Trac® ATX with dual rear strut suspension	Advanced Active-Trac® with independent rear castors
Maximum speed ¹	Up to 6 km/h (4 mph)	Up to 8 km/h (5 mph)	Up to 9.5 km/h (6 mph)	Up to 6 km/h (4 mph)
Brakes	Electronic, regenerative disc brakes	Electronic, regenerative disc brakes	Electronic, regenerative disc brakes	Electronic, regenerative disc brakes
Ground clearance	6.35 cm (2.5")	7.62 cm (3")	8.89 cm (3.5")	6.35 cm (2.5")
Turning radius	54 cm (21.25")	51 cm (20")	57 cm (22.5")	52 cm (20.5") (w/o foot riggings)
Overall size	Length: 88.27 cm (34.75") Width: 59 cm (23.25")	Length: 92.71 cm (36.5") Width: 58.42 cm (23")	Length: 91.44 cm (36") Width: 25.875" (65.72 cm)	Length: 93 cm (36.6") Width: 58.42 cm (23")
Drive train	Mid-Wheel 6® 2-pole Glide	Mid-Wheel 6® 2-pole Glide motor	Mid-Wheel 6® 4-pole 9.5 km/h (6 mph) motor	Mid-Wheel Drive 2-pole Glide
Batteries ¹	Two 12 volt, U-1	Two 12 volt, NF-22	Two 12 volt, NF-22, Group 24 or Group 34	Two 12 volt, U-1
Range per charge	Up to 24 km (15 miles)	Up to 32 km (20 miles)	Up to 32 km (20 miles)	Up to 24 km (15 miles)
Battery charger	5A, Off-board	8A, Off-board	8A, Off-board	5A, Off-board
Electronics	80A, Q-Logic	80A, Q-Logic	100A, Q-Logic	80A, Q-Logic
Weight capacity	136 kg (300 lbs.)	136 kg (300 lbs.)	136 kg (300 lbs.)	72.5 kg (160 lbs.)
Base weight	51 kg (112 lbs.)	58.5 kg (129 lbs.)	67 kg (147 lbs.)	51 kg (112 lbs.)
Battery weight ²	11 kg (24.5 lbs.) (each)	17 kg (38 lbs.) (each)	17 kg (38 lbs.) (each) or 24 kg (53.5 lbs.) (each)	11 kg (24.5 lbs.) (each)

Footnotes:

- 1) Range and speed vary with user weight, terrain type, battery charge, battery condition and tyre pressure.
- 2) Battery weight may vary based on manufacturer.

WARRANTY: 2-year limited warranty


Colours available on Q600 and Q610



Colours available on Q600, Q610 and Dynamo 1107



Colours available on Q6000Z



Colours available on Q600

