

# Pride® QUEST

FOLDING SCOOTER

THE **ULTIMATE**  
COMPACT, **FOLDABLE**  
TRAVEL **SCOOTER**

Ready for transport and easily folded, the 4-wheel **Quest Folding Scooter** is perfect for the active lifestyle. Offering great performance indoors and out at a maximum speed of 6 km/h (4 mph) and a range per charge of up to 20 km (12.5 miles), the Quest Folding Scooter provides a comfortable and convenient ride. It is super portable and is a great space saver that adjusts to fit perfectly inside any closet, corner, or vehicle trunk. The Quest Folding Scooter is sleek, stylish, and extremely easy to operate, making it the perfect travel mobility scooter choice.



Scooter Folded



Easy to Transport



PRIDE® MOBILITY SCOOTERS  
LIVE YOUR BEST®

[www.pridemobility.com.au](http://www.pridemobility.com.au)

# SPECIFICATIONS

## QUEST FOLDING SCOOTER

4W: S19

Weight capacity:	115 kg (250 lbs.)
Maximum speed <sup>1</sup> :	Up to 6 km/h (4 mph)
Ground clearance <sup>2</sup> :	6 cm (2.3")
Turning radius <sup>2</sup> :	82 cm (32")
Overall length <sup>2</sup> :	93 cm (37")
Overall width <sup>2</sup> :	48.5 cm (19")
Drive wheels <sup>2</sup> :	5 cm x 20 cm (2" x 8") solid
Front wheels <sup>2</sup> :	4 cm x 18 cm (1.6" x 7") solid
Range per charge <sup>1,4</sup> :	Up to 20 km (12.5 miles) (lithium-ion);
Total weight without batteries <sup>6</sup> :	23.5 kg (51.4 lbs.)
Heaviest piece when disassembled:	23.5 kg (51.4 lbs.) (unit does not disassemble - unit weight without battery pack)
Standard seating:	Type: Foldable, compact, vinyl-covered molded plastic Material: Gray Vinyl Dimensions: 41 cm x 33 cm (16" x 13")
Battery requirements <sup>5</sup> :	(1) 24-volt (lithium-ion) Size: (1) 11.5 AH (lithium-ion) Weight <sup>3</sup> : 11.5 AH lithium-ion- 2.3 kg (5 lbs.) each (1 required)
Battery charger:	2.5-amp, off-board

Colour:




White

## FEATURES

- Unique folding design
- Folds in just three easy steps
- Adjustable tiller
- Up to 20 km (12.5 miles) range per charge
- Easy to transport at home, hotels, or airports
- Can be stowed upright or laid down
- Articulated front end adds greater stability
- Front and rear LED lights

TO FOLD:



Fold Seat




Pull Up Lever to Fold




Fold and Lock Tiller

TO UNFOLD:



Pull Lever



Unfold Scooter



Unfold Seat, Lock in Place



Note: Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Go-Go® Folding Scooter and its lithium-ion batteries are permitted.

- <sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- <sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- <sup>3</sup> AGM or gel cell type required.
- <sup>4</sup> Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance test conducted at maximum weight capacity.
- <sup>5</sup> Battery weight may vary based on manufacturer.

All specifications provided by Heartway® Medical Products.  
Product distributed by Pride Mobility®.

**NOTE: All specifications subject to change without notice.**

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.



**Pride**  
Mobility Products  
Australia Pty. Ltd.  
— LIVE YOUR BEST —

PRIDE® MOBILITY SCOOTERS  
LIVE YOUR BEST®

(Tel) 03 8770 9600 • (Fax) 03 9703 2720

20-24 Apollo Drive Hallam, Victoria 3803

ACN 088 609 661 ABN 66 088 609 661

[www.pridemobility.com.au](http://www.pridemobility.com.au)

