

Unequaled Versatility

Fusion

Rear-Wheel Drive Power Wheelchair



A rugged design meets precise manoeuvrability

with the Fusion. This rear-wheel drive power chair delivers aggressive, terrain-gripping torque and comes standard with a top speed of up to 10 km/h (6 mph). The innovative, compact design of the Fusion makes it highly manoeuvrable in tight spaces without sacrificing traditional rear-wheel drive performance.

SPECIFICATIONS

Weight capacity	136 kg (300 lbs.)
Maximum speed¹	Up to 10 km/h (6 mph) (standard)
Ground clearance²	7.3 cm (3")
Turning radius²	63.5 cm (25")
Overall length²	89 cm (35")
Overall width²	63.5 cm (25") (at rear wheels)
Drive wheels	35.5 cm (14") pneumatic
Anti-tips	6.35 cm (2.5") solid (rear)
Casters	23 cm (9") pneumatic (front)
Suspension	N/A
Drive train	Two-motor, rear-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	70 amp PG VR2 controller
Standard seat	Sling (seat integrated with base)
Specialty controls	No
Battery charger	Off-board, 5 amp
Per-charge range¹	Up to 25 km (15.5 miles)
Battery requirements	(2) 12 volt deep-cycle
Battery size (largest)	50 AH
Battery weight³	17 kg (37 lbs.) each
Base weight	76.2 kg (170 lbs.) (w/o batteries)
Standard seat weight	17.23 kg (38 lbs.) (seat is integrated with base)
Std. seating dimensions	40.6 cm x 45.7 cm (16" x 18")
Warranty	2-year limited

COLOUR



Silver

Pride[®]
 Mobility Products
 Australia Pty. Ltd.

20-24 Apollo Drive Hallam, Victoria 3803
 ACN 088 609 661 ABN 66 088 609 661

(Tel) 03 8770 9600 • (Fax) 03 9703 2720 • www.pridemobility.com.au



The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

FEATURES

- Outstanding tight space manoeuvrability
- Maximum speed up to 10 km/h (6 mph)
- High-performance, high efficiency motor package
- Unique design seating system is integrated with the base
- 2-pole motors
- Power tilt

Fusion

shown in tilted position
 in Silver



FOOTNOTES:

- 1) Range and speed vary with user weight, terrain type, battery charge, battery condition and tyre pressure.
- 2) Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
- 3) Battery weight may vary based on manufacturer.

