


	Revo [®] 3-wheel	Revo [®] 4-wheel
Model number:	SC63	SC64
Weight capacity:	300 lbs.	300 lbs.
Maximum speed:	Up to 4.5 mph	Up to 4.5 mph
Ground clearance:	2.25"	2.25"
Turning radius:	38"	50"
Overall length:	43"	44.75"
Overall width:	21.5"	21.5"
Wheels		
Front:	2" x 8" solid	2" x 8" solid
Rear:	3" x 9" solid	3" x 9" solid
Anti-tips:	Rear	Rear
Suspension:	N/A	N/A
Transport weight		
Total w/o batteries:	106.5 lbs.	125 lbs.
Heaviest piece when disassembled:	48 lbs. (rear section)	51.5 lbs. (front section)
Drive train:	Sealed transaxle, 24 volt DC motor	
Dual braking system:	Electronic, regenerative and electromechanical (electronic, regenerative disc brakes)	
Power elevating seat:	N/A	
Std. Seating:	Compact foldable molded plastic	
Battery charger:	O -board, 5 amp (standard)	
Per-charge range:	Up to 20 miles	
Battery requirements:	(2) 12V deep cycle (AGM or gel cell type recommended)	
Battery size:	U-1	
Warranty:	3-year limited warranty	

Colors:



Red

Blue

Benefits


- Exclusive black non-scuffing tires
- Innovative interlocking drop-in battery packs for easy removal/installation
- Easily disassembles into 5 lightweight, easy-to-manage pieces
- 300 lbs. weight capacity
- State-of-the-art, microprocessor-based motor controller
- Seat offers height adjustment, 360° rotation, and flip-up armrests for optimal user comfort
- Adjustable tiller for maximum comfort
- Headlight and taillight for added safety
- Easily accessible charger port and circuit breaker

Options

- Forearm crutch holder
- Lap belt
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Weather cover
- Saddlebag
- Quad cane holder
- Cane, clip attachment

Pride SilverStar Lifts & Ramps 

- Backpacker[®] Plus
- Full Platform
- Silver Boom 300



Full platform



Backpacker[®] Plus

These lifts are compatible with the Revo[®].

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

Pride[®]
Mobility Products Corp.

REVO®

THE REVO PROVIDES SLEEK STYLING
AND SUPER PORTABILITY.

REVO®

The Revo® delivers full-sized scooter style, performance, and features in a compact, lightweight design that easily disassembles for transport while remaining highly affordable.



Pride®
Mobility Products Corp.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

REVOSS04-08