

Stylus[®]



Quality and Comfort...for Life.

The Stylus[®] is designed to meet the needs of active users. It's extremely lightweight, highly durable, and provides enhanced fit and function. The Stylus offers all this in a very cost effective package.

PRIDE[®]
PERFORMANCE MOBILITY
— LIVE YOUR BEST[™] —

SPECIFICATIONS

PRODUCT WEIGHT	33.5 lbs. without leg rests ¹
WEIGHT CAPACITY	250 lbs.
FRAME MATERIAL	Lightweight aluminum
SEAT-TO-FLOOR HEIGHT²	16"-20" heights available
SEAT WIDTHS	16", 18" and 20"
SEAT DEPTHS	16" and 18"
BACK HEIGHTS	17" - 19"
REAR WHEELS	22" or 24" composite mag wheels with polyurethane tires and aluminum handrims (standard)
WHEEL LOCKS	Push-to-lock
CASTERS	5", 6" and 8" casters available (flat-free)
WARRANTY	Lifetime limited on side-frames and cross-frames; 1-year limited on materials

FEATURES

- Adjustable tension seat upholstery increases comfort and positioning
- Multiple-position rear wheel and caster mounts optimize center-of-gravity and seat-to-floor height adjustments
- Quick-release, removable rear wheels
- Full-profile 22" or 24" x 1 3/8" rear wheels for enhanced outdoor capabilities
- Long handle wheel locks ensure comfort and security
- Industry-compatible leg rests with aluminum footplates for durability

OPTIONS

- Synergy[®] Cushions
- Synergy[®] Shape Back
- Angle-adjustable back canes
- Tension adjustable back upholstery
- Adjustable, removable rear anti-tips
- Oxygen holder
- Elevating leg rests with calf pads
- Height-adjustable armrests, desk or full length
- 6" x 2" flat-free casters
- 8" x 1 1/4" pneumatic casters
- Pull-to-lock wheel locks
- Wheel lock extensions
- 24" pneumatic rear wheels
- Vinyl-coated handrims
- Nylon backpack

RECOMMENDED LIFTS

Silver Star[®]
Lifts & Ramps

- Caddy
- Silver Boom 300



STYLUS
shown with Synergy[®] Solution[®] cushion

Pride[®]
Mobility Products Co.
— LIVE YOUR BEST[™] —

St. Catharines, Ontario ■ 888-570-1113
Visit us on the web: www.pridemobility.com

The information contained herein is correct at the time of publication.
We reserve the right to alter specifications without prior notice.

¹) Standard 16" x 16" seat size and height-adjustable, full length arms

²) Height varies from 16" - 20" based on wheel size